



# FavorPrep™ Soil RNA Isolation Mini Kit

Cat.No. : FASRK 004, 4 preps  
 FASRK 050, 50 Preps  
 FASRK 100, 100 Preps

**(For Research Use Only)**

## Kit Contents:

	FASRK 004 (4 preps_sample)	FASRK 050 (50 preps)	FASRK 100 (100 preps)
Glass Beads	4 vials	50 vials	100 vials
SR1 Buffer	5 ml	60 ml	120 ml
SR2 Buffer	0.3 ml	4 ml	8 ml
SR3 Buffer	0.6 ml	8 ml	16 ml
SR4 Buffer	1.5 ml	20 ml	40 ml
Wash Buffer (concentrate) *	1.5 ml	15 ml	30 ml
RNase-free water	1.5 ml	20 ml	40 ml
RNA Mini Column	4 pcs	50 pcs	100 pcs
Collection Tube	8 pcs	100 pcs	200 pcs
Elution Tube (1.5 ml tube)	4 pcs	50 pcs	100 pcs
2.0 ml tube	4 pcs	50 pcs	100 pcs
User Manual	1	1	1

* Preparation of Wash Buffer by adding ethanol (96% ~100%) for first use:			
Cat. No:	FASRK 004	FASRK 050	FASRK 100
ethanol volume for Wash Buffer	6 ml	60 ml	120 ml

## Specification:

Principle: spin column (silica membrane)  
 Sample: 0.25 ~ 1 g  
 Operation time: < 60 min  
 Elution volume: 40 µl


## Important Notes:

1. Make sure everything is RNase-free when handling RNA.
2. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
3. Add indicated volume of ethanol (96-100%) to Wash Buffer before use.
4. Phenol and chloroform are hazardous to human health, perform the procedures involving phenol or chloroform in a chemical fume hood.
5. Preheat RNase free water to 60 °C for elution step.

## Additional material to be provided by user:

- Ethanol (RNase-free, 96~100%)
- Isopropanol
- Water-saturated phenol
- Chloroform
- Microcentrifuge for 1.5 ml and 2.0 ml tube capable of at least 13,000 x g
- Centrifuge for 15 ml tube capable of at least 4,000 X g
- Water Bath or Dry Bath
- 1.5 ml centrifuge tube
- 15 ml centrifuge tube

## Brief Procedure:



## General Protocol:

### Please Read Important Notes Before Starting Following Steps.

1. Add the Glass Beads (provided) and 1 g of soil sample to a 15 ml centrifuge tube (not provided).
  - If the sample is liquid, add 500  $\mu$ l of sample to the 15 ml Beads Tube.
2. Add 1 ml of SR1 Buffer to the sample, vortex at maximum speed for 1 minute.
3. Add 50  $\mu$ l of SR2 Buffer to the sample, vortex at maximum speed for 1 minute.
4. Add 1 ml of water saturated-phenol and 200  $\mu$ l of chloroform to the sample, mix well by plus-vortexing for 5 minutes to lysis the sample.
  - Caution!** Phenol and chloroform are hazardous, performing Step 4 and Step 5 in a chemical fume hood.
5. Centrifuge at 4,000 x g for 10 minutes to form aqueous and organic phase.
6. Carefully transfer the upper aqueous phase to a 2.0 ml centrifuge tube (provided). Measure the volume of the upper aqueous.
  - Carefully transfer the upper aqueous phase and do not disrupt the inter phase that consist of phenol and chloroform.
  - Avoid pipetting any debris and pellet.
7. Add 1 volume of isopropanol, vortex to mix well. centrifuge at full speed (13,000 x g) for 10 min to form a pellet.
8. Carefully discard the supernatant and invert the tube on the paper towel for 5 min to remove residual liquid.
  - Do not disrupt the pellet when discard the supernatant.
  - Depending on the soil type, the dark color of pellet is not consist .
9. Add 200  $\mu$ l of RNase-free ddH<sub>2</sub>O, vortex to dissolve the pellet completely.
  - Incubate the tube at 45°C for 10 min if the pellet is hard to dissolve.
10. Add 100  $\mu$ l of SR3 Buffer to the sample, mix well by vortexing. Incubate the sample at room temperature for 2 minutes.
  - Note: SR3 Buffer must be suspended completely by vigorously vortexing before every using.
  - use 1ml pipettor and cut off the end of 1 ml tip to make it easier for pipetting the SR3 Buffer.
11. Centrifuge at full speed (13,000 x g) for 2 minutes.
12. Carefully transfer the clarified lysate to a 1.5 ml microcentrifuge (not provided).
  - Avoid pipetting any debris and pellet.
13. Add 250  $\mu$ l of SR4 Buffer and 250  $\mu$ l of ethanol (RNase-free, 96~100%) to the clarified lysate, mix thoroughly by pulse-vortexing.
14. Place a RNA Mini Column into a Collection Tube and transfer all of the sample mixture to the RNA Column. Centrifuge at full speed (13,000 x g) for 1 min then discard the flow-through. Place the RNA Column to a new Collection Tube.
15. Add 650  $\mu$ l of Wash Buffer (ethanol added) to RNA Mini Column. Centrifuge at full speed (13,000 x g) for 1 min then discard the flow-through.
  - Make sure that ethanol (RNase-free, 96~100%) has been added into Wash Buffer when first open.
16. Repeat step 15 for one more time.
17. Centrifuge at full speed (13,000 x g) for an additional 3 min to dry the RNA Mini Column.
  - Important step! This step will avoid the residual liquid to inhibit subsequent enzymatic reactions.
18. Place RNA Column into a Elution Tube, Add 40  $\mu$ l of pre-heated RNase-free ddH<sub>2</sub>O to the membrane center of the RNA Mini Column. Stand the RNA Column for 2 min at room temperature.
  - Important step! For effective elution, make sure that the RNase-free ddH<sub>2</sub>O is dispensed onto the membrane center and is absorbed completely.
19. Centrifuge at full speed (13,000 x g) for 1 min to elute RNA.