
SUB Series Horizontal Electrophoresis Units

SUB6, SUB10, SUB13, SUB15, SUB20, SUB25,
SUB20C, SUB25C and SUBHT

Page Finder

Important Information	ii
Waste Electrical & Electronic Equipment (WEEE)	ix
Specifications.....	1
Safety Precautions.....	2
General Care and Maintenance	3
Gel Loading Instructions for the SUBHT	5
Cleaning and Storage.....	6
Packing Lists	7
Operational	9
SUB6 Mini Horizontal	9
SUB10 Mini-Plus Horizontal	10
SUB13 Midi Horizontal	11
SUB15 Standard Horizontal	12
SUB20 Maxi-Standard Horizontal.....	13
SUB25 Maxi-Plus Standard Horizontal	14
SUB20C Maxi-Cooled Horizontal.....	15
SUB25C Maxi-Plus Cooled Horizontal.....	16
SUBHT High Throughput Horizontal.....	17
Using The Horizontal Gel Electrophoresis Units.....	18
Preparation of Agarose Gels	18
Agarose Gel Electrophoresis.....	19
Gel Pouring.....	20
Running the Gel	22

Important Information – English

- If this equipment is used in a manner not specified by Hoefer, Inc. the protection provided by the equipment may be impaired.
- This instrument is designed for indoor laboratory use only.
- Only accessories and parts approved or supplied by Hoefer, Inc. may be used for operating, maintaining, and servicing this product.
- Only use a power supply that is CE marked or safety certified by a nationally recognized testing laboratory.
- The safety lid must be in place before connecting the power supply leads to a power supply.
- Turn all power supply controls off and disconnect the power leads before removing the safety lid.
- Circulate only water or 50/50 water/ethylene glycol through the heat exchanger if so equipped. Do not connect the heat exchanger to a water tap or any coolant source where the water pressure is unregulated.
- Never introduce antifreeze or any organic solvent into any part of the instrument. Organic solvents will cause irreparable damage to the unit!
- Do not operate with buffer temperatures above the maximum specified technical specifications. Overheating will cause irreparable damage to the unit!

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Waste Electrical and Electronic Equipment (WEEE)

English

This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.

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Specifications

Application

For laboratory use only in the electrophoretic separation of proteins and nucleic acids.

Country of Origin: United States

Manufacturer:

Hoefer, Inc.
84 October Hill Road
Holliston, MA 01746
Toll Free: 1-800-227-4750
Phone: 1-508-893-8999
Fax: 1-508-429-5732
E-mail: sales@hoeferinc.com

Construction

- Rugged acrylic construction.
- All acrylic joints chemically bonded.
- Doubly insulated cables, rated safe up to 1,000 V.
- Gold plated electrical connectors, corrosion-free and rated safe up to 1,000 V.
- Recessed power connectors, integral with the safety lid.
- 0.2 mm diameter platinum electrodes, 99.99% pure.
- User replaceable platinum electrodes.
- User-friendly clamping system.
- Wide range of accessories.

Environmental Conditions

- This apparatus is intended for indoor use only.
- This apparatus can be operated safely up to an altitude of 2,000 m.
- The normal operating temperature range is between 4°C and 65°C.
- Maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.
- The apparatus is rated POLLUTION DEGREE 2 in accordance with IEC 664. POLLUTION DEGREE 2, states that: "Normally only non-conductive pollution occurs. Occasionally, however, a temporary conductivity caused by condensation must be expected".

Safety Precautions

- READ the instructions before using the apparatus.
- The exclamation mark within an equilateral triangle is intended to alert the user to the presence of important operating and maintenance instructions in the literature accompanying the instrument.
- Always unplug electrophoresis units from their power supply before removing the safety cover. Unplug the power supply from the mains **FIRST**, then disconnect the leads.
- **DO NOT** exceed the maximum operating voltage or current.
- **DO NOT** operate the electrophoresis units in metal trays.
- It is advisable to always place the unit on a plastic safety tray during electrophoresis.
- Ethidium bromide is a known mutagen. Wear effective protective clothing when handling and follow recommended handling and disposal procedures.
- Ultraviolet light can cause permanent damage to the retina and skin cancer. Wear UV safety goggles and cover any exposed skin when using a UV light source.
- Following the replacement of a platinum electrode, have the unit inspected and approved by your safety officer prior to use.
- **DO NOT** fill the unit with running buffer above the maximum fill lines.
- **DO NOT** move the unit when it is running.
- **CAUTION:** During electrophoresis, very low quantities of various gases are produced at the electrodes. The type of gas produced depends on the composition of the buffer employed. To disperse these gases make sure that the apparatus is run in a well ventilated area.

IMPORTANT! Acrylic plastic is NOT resistant to aromatic or halogenated hydrocarbons, ketones, esters, alcohol's (over 25%) and acids (over 25%), they will cause "crazing" especially of the UV transparent plastic and should NOT be used for cleaning. DO NOT use abrasive creams or scourers. Dry components with clean laboratory tissues prior to use.

General Care and Maintenance

- To remove the safety lid, push thumbs down on the plastic lugs and lift the lid vertically with your fingers.
- Before use, clean and dry the apparatus with DISTILLED WATER ONLY.
- Before use, and then on a monthly basis, check the unit for any leaks at the bonded joints. Place the unit on a sheet of dry tissue and then fill with DISTILLED WATER ONLY to the maximum fill line. Any leakage will be seen on the tissue paper. If any leakage is seen DO NOT ATTEMPT TO REPAIR OR USE THE APPARATUS, but notify Hoefor immediately.
- The replacement platinum electrodes are partially shrouded for protection. However, when cleaning the main tank DO NOT use cleaning brushes in the electrode area. Usually a thorough rinse with distilled water is all that is required.
- Ensure that the connectors are clean and dry before usage or storage.

Storage of Water Cooled Units

We recommend complete drainage of cooling plate coil units before storage. Store in a dark cupboard or cold room.

Using the Cooling Platform

The Cooling Platform will already contain a small quantity of water from control tests. The Cooling Platform can be used in two ways. Static water can be used as a simple heat sink or the tank can be actively regulated using flowing water from a tap or water bath/flow cooler.

Static Temperature Regulation


- 1 Attach a short length of rubber hose to each connector.
- 2 Incline the unit at an angle of approximately 45 degrees with the ports uppermost.
- 3 Use a funnel to fill the Cooling Platform with deionized water containing 0.02% (w/v) sodium azide (preservative to prevent bacterial and algal growth).
- 4 When filled, keep the unit inclined and attach clamps to each piece of rubber hose.
- 5 The unit can be cooled before an electrophoresis run if required. **DO NOT FREEZE.**

Active Temperature Regulation

- 1 Attach a short length of rubber hose to each connector.
- 2 Attach one end of the rubber hose to the outlet port of a circulating water bath/flow cooler and the other end to the inlet port. Alternatively attach one end of a rubber hose to a water supply and allow the other rubber hose drain to waste.
- 3 **The maximum recommended water flow rate is 1 Liter/min. DO NOT exceed this figure.**
- 4 If you are using a circulating water bath/flow cooler, which exceeds this flow rate, you can attach a T-connector in line. One branch of the connector can return water to the bath and the other can flow to the Cooling Platform and incorporate a flow regulator such as an adjustable tubing clamp. Measure and adjust the flow rate before attaching the line to the gel unit.

Gel Loading Instructions for the SUBHT

- 1 Use a 12 channel multi-pipette to aspirate sample from row A of a 96-well plate (Wells A1 – A12).
- 2 Place the left-most tip (containing sample A1) into the corresponding gel well labeled A1 (4th well in from the left side) and dispense. Samples from wells A1 though A12 will be dispensed in every other well of the gel.
- 3 Aspirate samples from row B (wells B1 – B12) and place left-most tip (containing sample B1) into the corresponding gel well labeled B1 (5th well in from the left side, to the right of well A1).
- 4 Repeat steps 1–3 for rows C and D, E and F, and G and H.
- 5 Use a single channel pipette to load up to 3 size marker samples on either side of each row.

IMPORTANT! Acrylic plastic is NOT resistant to aromatic or halogenated hydrocarbons, ketones, esters, alcohol's (over 25%) and acids (over 25%), they will cause "crazing" especially of the UV transparent plastic and should NOT be used for cleaning. DO NOT use abrasive creams or scourers. Dry components with clean tissues prior to use.

At the End of the Run

- 1 Turn the power supply settings to zero, turn off mains supply and disconnect the power leads.
- 2 Visualize the run progression or final separation on a UV transilluminator.
- 3 After each use, rinse the apparatus with **DISTILLED WATER ONLY.**
- 4 After each use, wash all components with a suitable water-soluble disinfectant, to reduce the risk of infection.
- 5 Ensure that the connectors are clean and dry before usage or storage.

Cleaning and Storage

After use, thoroughly rinse all components in deionized water and gently dry the gold electrode connectors with a soft tissue.

NEVER USE ORGANIC SOLVENTS.

Packing Lists

SUB6

No. Items	Description	Catalog No.
1	Running tank + lid	—
1	Gel casting tray	SUB6-UT
2	Combs, 1.0 mm 8 sample	SUB6-C1-8
6	Colored loading strips	SUB6-CS

SUB10

1	Running tank + lid	—
1	Gel casting tray	SUB10-UT
2	Combs, 1.0 mm 16 sample	SUB10-C1-16
6	Colored loading strips	SUB10-CS

SUB13

1	Running tank + lid	—
1	Gel casting tray	SUB13-UT
2	Combs, 1.0 mm 16 sample	SUB13-C1-16
2	Buffer recirculation ports	SUB-BRP
6	Colored loading strips	SUB13-CS

SUB15

1	Running tank + lid	—
1	Gel casting tray	SUB15-UT
2	Casting gates + silicone seals	SUB15-CG
2	Combs, 1.0 mm 16 sample	SUB15-C1-16
2	Buffer recirculation Ports	SUB-BRP
6	Colored loading strips	SUB15-CS

SUB20

1	Running tank + lid	—
1	Gel casting tray	SUB20-UT
2	Casting gates + silicone seals	SUB20-CG
2	Combs, 1.0 mm 16 sample	SUB20-C1-16
2	Buffer recirculation ports	SUB-BRP
6	Colored loading strips	SUB20-CS

SUB25

No. Items	Description	Catalog No.
1	Running tank + lid	—
1	Gel casting tray	SUB25-UT
2	Casting gates + silicone seals	SUB25-CG
6	Combs, 1.0 mm 26 sample	SUB25-C1-26MC
2	Buffer recirculation ports	SUB-BRP
6	Colored loading strips	SUB25-CS

SUB20C

1	Running tank + lid	—
1	Gel casting tray	SUB20-UT
2	Casting gates + silicone seals	SUB20-CG
2	Combs, 1.0 mm 16 sample	SUB20-C1-16
2	Buffer recirculation ports	SUB-BRP
6	Colored loading strips	SUB20-CS


SUB25C

1	Running tank + lid	—
1	Gel casting tray	SUB25-UT
2	Casting gates + silicone seals	SUB25-CG
6	Combs, 1.0 mm 26 sample	SUB25-C1-26MC
2	Buffer recirculation ports	SUB-BRP
6	Colored loading strips	SUB25-CS

SUBHT

1	Running tank + lid	—
1	Labeled gel casting tray	SUBHT-LUT
2	Casting gates + silicone seals	SUBHT-CG
2	Combs, 1.0 mm 30 sample	SUB15-C1-30MC
2	Buffer recirculation ports	SUB-BRP
6	Colored loading strips	SUBHT-CS

Operational

SUB6 Mini Horizontal

The SUB6 Mini Horizontal gel electrophoresis unit is ideal for routine preparatory and analytical electrophoresis techniques.


Technical Specifications

Unit dimensions (W × L × H)	13 × 24 × 6.5 cm
Gel dimensions (W × L)	6 × 7.5 cm
Buffer volume	325 ml
Maximum sample capacity	32
Combs	2
Comb thickness	1, 1.5 or 2 mm
Comb throughput	8 to 16 samples
Comb slots	2
Migration distance between comb slots	3.5 cm
Recommended running voltage	70 to 90 V

SUB6 Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (µl)
SUB6-C1-8	1	8	4.5	2	20
SUB6-C1-12MC	1	12*	2.5	2	11
SUB6-C1-16	1	16	2.2	1	10
SUB6-C1.5-8	1.5	8	4.5	2	30
SUB6-C1.5-12MC	1.5	12*	2.5	2	17
SUB6-C1.5-16	1.5	16	2.2	1	15
SUB6-C2-8	2	8	4.5	2	40
SUB6-C2-12MC	2	12*	2.5	2	22
SUB6-C2-16	2	16	2.2	1	20

*Multi-channel pipette compatible

SUB10 Mini-Plus Horizontal

The SUB10 Mini-Plus Horizontal gel electrophoresis unit is ideal for routine preparatory and analytical electrophoresis techniques.


Technical Specifications

Unit dimensions (W × L × H)	16.5 × 23 × 6.5 cm
Gel dimensions (W × L)	10 × 11.5 cm
Buffer volume	450 ml
Maximum sample capacity	80
Combs	2
Comb thickness	1, 1.5 or 2 mm
Comb throughput	8 to 20 samples
Comb slots	4
Migration distance between comb slots	2.5 cm
Recommended running voltage	75 to 125 V

SUB10 Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (µl)
SUB10-C1-10MC	1	10*	7	2	30
SUB10-C1-12	1	12	5.5	2	25
SUB10-C1-16	1	16	3.6	2	15
SUB10-C1-20MC	1	20*	3	1.5	12
SUB10-C1.5-8	1.5	8	9	2.5	60
SUB10-C1.5-10MC	1.5	10*	7	2	45
SUB10-C1.5-12	1.5	12	5.5	2	35
SUB10-C1.5-16	1.5	16	3.6	2	25
SUB10-C1.5-20MC	1.5	20*	3	1.5	20
SUB10-C2-8	2	8	9	2.5	80
SUB10-C2-10MC	2	10*	7	2	60
SUB10-C2-12	2	12	5.5	2	50
SUB10-C2-16	2	16	3.6	2	30
SUB10-C2-20MC	2	20*	3	1.5	25

*Multi-channel pipette compatible

SUB13 Midi Horizontal

The SUB13 Midi Horizontal gel electrophoresis unit is ideal for analytical and preparative studies of nucleic acids.


Technical Specifications

Unit dimensions (W × L × H)	20 × 32 × 7 cm
Gel dimensions (W × L)	12.8 × 15 cm
Buffer volume	900 ml
Buffer re-circulation ports	2
Maximum sample capacity	112
Combs	2
Comb thickness	1, 1.5 or 2 mm
Comb throughput	10 to 28 samples
Comb slots	4
Migration distance between comb slots	3.5 cm
Recommended running voltage	100 to 125 V

SUB13 Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (µl)
SUB13-C1-10	1	10	9.5	2.5	40
SUB13-C1-12MC	1	12	8	2	35
SUB13-C1-16	1	16	5.5	2	25
SUB13-C1-20	1	20	4	2	17
SUB13-C1-24	1	24	3	2	13
SUB13-C1-28MC	1	28*	3	1.5	13
SUB13-C1.5-10	1.5	10	9.5	2.5	60
SUB13-C1.5-12MC	1.5	12	8	2	50
SUB13-C1.5-16	1.5	16	5.5	2	35
SUB13-C1.5-20	1.5	20	4	2	25
SUB13-C1.5-24	1.5	24	3	2	20
SUB13-C1.5-28MC	1.5	28*	3	1.5	20
SUB13-C2-10	2	10	9.5	2.5	85
SUB13-C2-12MC	2	12	8	2	70
SUB13-C2-16	2	16	5.5	2	50
SUB13-C2-20	2	20	4	2	35
SUB13-C2-24	2	24	3	2	25
SUB13-C2-28MC	2	28*	3	1.5	25

*Multi-channel pipette compatible

SUB15 Standard Horizontal

The SUB15 Standard Horizontal gel electrophoresis unit provides the comb-throughput and tray-size options most popular in today's market.


Technical Specifications

Unit dimensions (W × L × H)	21.5 × 33.5 × 7 cm
Gel size (W × L)	15 × 15 cm
Buffer volume	1200 ml
Buffer re-circulation ports	2
Maximum sample capacity	120
Combs	2
Comb thickness	1, 1.5 or 2 mm
Comb throughput	1 to 30 samples
Comb slots	4
Migration distance between comb slots	3.5 cm
Recommended running voltage	100 to 125 V

SUB15 Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (μl)
SUB15-C1-10	1	10	12.8	2	64
SUB15-C1-16MC	1	16*	6.4	2	32
SUB15-C1-20	1	20	5.9	2	29
SUB15-C1-25	1	25	3.5	2	17
SUB15-C1-30MC	1	30*	2.5	2	12
SUB15-C1.5-10	1.5	10	12.8	2	96
SUB15-C1.5-16MC	1.5	16*	6.4	2	48
SUB15-C1.5-20	1.5	20	5.9	2	43
SUB15-C1.5-25	1.5	25	3.5	2	25
SUB15-C1.5-30MC	1.5	30*	2.5	2	18
SUB15-C2-10	2	10	12.8	2	128
SUB15 -C2-16MC	2	16*	6.4	2	64
SUB15-C2-20	2	20	5.9	2	58
SUB15-C2-25	2	25	3.5	2	34
SUB15-C2-30MC	2	30*	2.5	2	24

*Multi-channel pipette compatible

SUB20 Maxi-Standard Horizontal

The SUB20 Maxi-Standard Horizontal gel electrophoresis unit is ideal for high-resolution analytical and preparative studies of nucleic acids.


Technical Specifications

Unit dimensions (W × L × H)	27 × 47.5 × 8 cm
Gel size (W × L)	20 × 20 cm
Buffer volume	2200 ml
Buffer re-circulation ports	2
Maximum sample capacity	160
Combs	2
Comb thickness	1, 1.5 or 2 mm
Comb throughput	16 to 40 samples
Comb slots	4
Migration distance between comb slots	5 cm
Recommended running voltage	150 to 175 V

SUB20 Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (µl)
SUB20-C1-16	1	16	8.5	3	35
SUB20-C1-20MC	1	20*	7	2	30
SUB20-C1-28	1	28	5	2	20
SUB20-C1-40MC	1	40*	3	2	13
SUB20-C1.5-16	1.5	16	8.5	3	55
SUB20-C1.5-20MC	1.5	20*	7	2	45
SUB20-C1.5-28	1.5	28	5	2	30
SUB20-C1.5-40MC	1.5	40*	3	2	19
SUB20-C2-16	2	16	8.5	3	75
SUB20-C2-20MC	2	20*	7	2	60
SUB20-C2-28	2	28	5	2	40
SUB20-C2-40MC	2	40*	3	2	25

*Multi-channel pipette compatible

SUB25 Maxi-Plus Standard Horizontal

The SUB25 Maxi-Plus Standard Horizontal gel electrophoresis unit is ideal for high throughput analysis of nucleic acids.


Technical Specifications

Unit dimensions (W × L × H)	33 × 56 × 9 cm
Gel size (W × L)	25 × 30 cm
Buffer volume	3000 ml
Buffer re-circulation ports	2
Maximum sample capacity	624
Combs	6
Comb thickness	1, 1.5 or 2 mm
Comb throughput	26 to 52 samples
Comb slots	12
Migration distance between comb slots	2 cm
Recommended running voltage	150 to 200 V

SUB25 Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (µl)
SUB25-C1-26MC	1	26*	7	2	30
SUB25-C1-52MC	1	52*	3	1.5	13
SUB25-C1.5-26MC	1.5	26*	7	2	45
SUB25-C1.5-52MC	1.5	52*	3	1.5	20
SUB25-C2-26MC	2	26*	7	2	60
SUB25-C2-52MC	2	52*	3	1.5	25

*Multi-channel pipette compatible

SUB20C Maxi-Cooled Horizontal

The SUB20C Maxi-Cooled Horizontal gel electrophoresis unit allows analytical and preparative studies of nucleic acids to be performed quickly at high voltage without compromising resolution.


Technical Specifications

Unit dimensions (W × L × H)	27 × 47.5 × 8 cm
Gel size (W × L)	20 × 20 cm
Buffer volume	2200 ml
Buffer re-circulation ports	2
Maximum sample capacity	160
Combs	2
Comb thickness	1, 1.5 or 2 mm
Comb throughput	16 to 40 samples
Comb slots	6
Cooling ports	2
Migration distance between comb slots	5 cm
Recommended running voltage	150 to 225 V

SUB20C Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (µl)
SUB20-C1-16	1	16	8.5	3	35
SUB20-C1-20MC	1	20*	7	2	30
SUB20-C1-28	1	28	5	2	20
SUB20-C1-40MC	1	40*	3	2	13
SUB20-C1.5-16	1.5	16	8.5	3	55
SUB20-C1.5-20MC	1.5	20*	7	2	45
SUB20-C1.5-28	1.5	28	5	2	30
SUB20-C1.5-40MC	1.5	40*	3	2	19
SUB20-C2-16	2	16	8.5	3	75
SUB20-C2-20MC	2	20*	7	2	60
SUB20-C2-28	2	28	5	2	40
SUB20-C2-40MC	2	40*	3	2	25

*Multi-channel pipette compatible

SUB25C Maxi-Plus Cooled Horizontal

The SUB25C Maxi-Plus Cooled Horizontal gel electrophoresis unit allows high throughput screening of nucleic acids to be performed quickly at high voltage without compromising resolution.


Technical Specifications

Unit dimensions (W × L × H)	33 × 56 × 9 cm
Gel size (W × L)	25 × 30 cm
Buffer volume	3000 ml
Buffer re-circulation ports	2
Maximum sample capacity	624
Combs	6
Comb thickness	1, 1.5 or 2 mm
Comb throughput	26 to 52 samples
Comb slots	12
Cooling ports	2
Migration distance between comb slots	2 cm
Recommended running voltage	150 to 250 V

SUB25C Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (µl)
SUB25-C1-26MC	1	26*	7	2	30
SUB25-C1-52MC	1	52*	3	1.5	13
SUB25-C1.5-26MC	1.5	26*	7	2	45
SUB25-C1.5-52MC	1.5	52*	3	1.5	20
SUB25-C2-26MC	2	26*	7	2	60
SUB25-C2-52MC	2	52*	3	1.5	25

*Multi-channel pipette compatible

SUBHT High Throughput Horizontal

The SUBHT High Throughput Horizontal gel electrophoresis unit is ideal for high throughput screening applications following PCR like Single Strand Polymorphism (SSP) analysis.

Technical Specifications

Unit dimensions (W × L × H)	21.5 × 33.5 × 7 cm
Gel size (W × L)	15 × 15 cm
Buffer volume	1200 ml
Buffer re-circulation ports	2
Maximum sample capacity	120
Combs	4
Comb thickness	1, 1.5 or 2 mm
Comb throughput	1 to 30 samples
Comb slots	4
Migration distance between comb slots	3.5 cm
Recommended running voltage	100 to 125 V

SUBHT Compatible Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (µl)
SUB15-C1-10	1	10	12.8	2	64
SUB15-C1-16MC	1	16*	6.4	2	32
SUB15-C1-20	1	20	5.9	2	29
SUB15-C1-25	1	25	3.5	2	17
SUB15-C1-30MC	1	30*	2.5	2	12
SUB15-C1.5-10	1.5	10	12.8	2	96
SUB15-C1.5-16MC	1.5	16*	6.4	2	48
SUB15-C1.5-20	1.5	20	5.9	2	43
SUB15-C1.5-25	1.5	25	3.5	2	25
SUB15-C1.5-30MC	1.5	30*	2.5	2	18
SUB15-C2-10	2	10	12.8	2	128
SUB15-C2-16MC	2	16*	6.4	2	64
SUB15-C2-20	2	20	5.9	2	58
SUB15-C2-25	2	25	3.5	2	34
SUB15-C2-30MC	2	30*	2.5	2	24

*Multi-channel pipette compatible

Using The Horizontal Gel Electrophoresis Units

Preparation of Agarose Gels

The amount of agarose needed to give a particular % gel needs to be calculated; for example a 0.8% gel can be made by dissolving 0.8 grams of agarose powder in 100 mls of running buffer. The type of running buffer used should always be the same for the gel as that used in the buffer tank. The agarose needs to be fully dissolved before the gel can properly form. This can be achieved by heating with swirling in a water bath or incubator set to 70 °C or by heating on a magnetic heating block with a magnetic stirring bar inserted. The flask should always be covered to prevent evaporation of the buffer and a higher concentration gel resulting. These methods may take longer than 1 hour for the agarose to fully dissolve. Alternatively, the agarose can be dissolved in about 5–10 minutes using a microwave oven. The agarose solution should be covered and the microwave set to low. The agarose dissolves better if the microwave is periodically stopped and the solution swirled. Before pouring, the agarose solution should be checked for undissolved agarose crystals which can affect the mobility characteristics of gels. If these are present, agarose dissolving should be continued. The agarose solution should be cooled to 50–60 °C before pouring. If poured at too hot a temperature, the gel will be more likely to leak and also the comb may become distorted.

Agarose Gel Electrophoresis

DNA Mobility

DNA fragments as small as 1 kb or less can be separated using agarose gel electrophoresis. For fragments smaller than 0.1 kb, polyacrylamide gels are more suited.

RNA Mobility

Either before or during electrophoresis, RNA should be denatured. For example,

- RNA fragments which have denatured with glyoxal and dimethyl sulphoxide can be separated on neutral agarose gels, or
- RNA can be fractionated on agarose gels containing methylmercuric hydroxide or formaldehyde.

RNA samples usually require longer runs or buffers that are easily depleted, so it is necessary to circulate the buffer. Northern analyses should not normally be run on a mini gel tank.

Separation Performance

Gel concentration, running buffer, voltage, temperature, conformation, and the presence of ethidium bromide all affect separation results.

Gel Concentration Selection

The range of fragment sizes to be separated will determine the choice of agarose concentration for a gel. Typical agarose concentration is 0.5% to 3.0%. For large DNA fragments low-percentage gels are required, while for small DNA fragments, high-percentage gels are recommended. Weak gels (<0.5% agarose) should be electrophoresed at low temperatures (e.g. ~4 °C). Agarose gels of 0.75% to 1.0%, for routine electrophoresis, are recommended for a wide range of separations (0.15 to 15 kb). 2–4% agarose gels are usually selected for PCR fragment resolution.

If the gel has to be subsequently photographed, thin gels (2–3 mm) with low-percentage agarose are better than thick or high-percentage gels. The latter produce increased opacity and autofluorescence.

Gel Pouring

SUB6, SUB10 and SUB13

Units with removable gel casting trays and casting within the tank.


IMPORTANT! Ensure that the agarose has cooled to between 50 and 60 °C to prevent apparatus distortion.

- 1 Make sure that the silicone rubber gasket is in position, i.e. pushed evenly into the groove at the ends of the gel casting tray.
- 2 Carefully place the casting tray sideways into the gel unit so that the gaskets form a seal with the sides of the running unit.
- 3 Push the casting tray down onto the gel platform.
- 4 If the casting tray is too tight, remove the silicone gaskets and replace it in the end grooves stretching it a little.
- 5 If the casting tray is too loose, remove the silicone gaskets and replace them in the end grooves compressing it a little.
- 6 Place the gel unit on a level surface.
- 7 Position the required comb(s) into the slot(s) in the gel casting tray.
- 8 Prepare 22.5, 60 or 100 mls of the desired % agarose for units SUB6, SUB10 or SUB13 respectively. This will give a gel depth of 5 mm. For a 10 mm deep gel, prepare 45, 125 or 200 mls respectively.
- 9 Using a 1 ml disposable pipette, seal the silicone seal with some of the cooled agarose.
- 10 Pour in the agarose smoothly so as not to generate bubbles. If bubbles do occur, these can be smoothed to the side of the gel and dispersed using a clean gloved hand.
- 11 Allow the agarose to set, ensuring that the gel remains undisturbed.
- 12 Carefully remove the gel casting tray from the running unit. Either remove the end silicone rubber gaskets, or leave them in place.
- 13 Place the gel casting tray in the running position such that the wells are nearest to the black electrode (cathode).
- 14 Fill the tank with running buffer so that the buffer just submerges the gel.
- 15 Carefully remove the comb.

SUB15, SUB20, SUB20C, SUB25, SUB25C and SUBHT

*Units with removable gel casting trays and
gel casting gates.*

IMPORTANT! Ensure that the agarose has cooled to between 50 and 60°C to prevent apparatus distortion.

- 1 Place the casting gates into the slots provided in the gel casting tray with the end silicone rubber gaskets facing outwards.
- 2 Place the gel unit on a level surface.
- 3 Position the required combs into the slots provided in the gel casting tray.
- 4 Prepare 115, 200 or 375 mls of the desired % agarose for units SUB15/SUBHT, SUB20/SUB20C, or SUB25/SUB25C respectively. This will give a gel depth of 5 mm. For a 10 mm deep gel, prepare 230, 400 or 750 mls respectively.
- 5 Using a 1 ml disposable pipette, seal the inside of the casting gates with some of the cooled agarose.
- 6 Pour in the agarose smoothly so as not to generate bubbles. If bubbles do occur, these can be smoothed to the side of the gel and dispersed using a clean gloved hand. Allow the agarose to set, ensuring that the gel remains undisturbed.
- 7 Carefully remove the gel casting gates and comb(s).
- 8 Place the gel casting tray in the running position such that the wells are nearest to the black electrode (cathode).
- 9 Fill the tank with running buffer so that the buffer just submerges the gel.

IMPORTANT! Do not exceed the recommended voltage or current as this may result in poor band resolution and may result in damage to the unit.

IMPORTANT! When re-circulating buffer, remember that the buffer flowing through the tubing is live. Take all necessary precautions. Warn other workers in the vicinity of the potential hazard. Seek the advice of your safety officer.

Running the Gel

- 1 Samples should be mixed with a suitable marker dye or loading buffer prior to loading. This aids sinking of the sample into the well, and visualization of how far the samples have run. Consult your laboratory manual for details on loading buffers.
- 2 Load the samples into the wells taking care not to damage the sides or bottoms of the wells. Replace the lid correctly **BEFORE** connecting the leads to the power supply.
- 3 Set the voltage and current to suit the electrophoretic application. As a guide, to obtain the optimum resolution of DNA fragments, agarose gels should not be run greater than 5V/cm.
- 4 Long runs may require buffer re-circulation, to prevent over heating and or buffer depletion. Re-circulation ports are provided in the SUB13, SUB15, SUB20, SUB20C, SUB25, SUB25C and SUBHT units.
- 5 The Cooled gel units have a built-in cooling block in the base. The side connectors can be attached, via tubing, to the mains water supply or to a circulating water bath.

Hoefler, Inc.

84 October Hill Road
Holliston, MA 01746

Toll Free: 1-800-227-4750

Phone: 1-508-893-8999

Fax: 1-508-893-0176

E-mail: support@hoeferinc.com

Web: www.hoeferinc.com

Hoefler is a registered trademark
of Hoefler, Inc.

© 2012 Hoefler, Inc. —
All rights reserved.

Printed in the USA.