


GENAXY CENTRIFUGE TUBE SAMPLE PESTLE

Centrifuge Tube Sample Pestle



- Material SP-T: Teflon, (Handle) Stainless Steel
- Material SP-P: Polypropylene
- Used to grind DNA pellets and/or cells in 1.5mL Microtubes
- Autoclavable

Cat. Number	Description	Qty. / Pack
GEN-SP-T	Centrifuge Tube Sample Pestle for 1.5mL Microtube, Teflon, 143mm long.	10
GEN-SP-P	Centrifuge Tube Sample Pestle for 1.5mL Microtube, Polypropylene, 85mm long.	100

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com