

# GENAXY HORIZONTAL SYSTEM

## MAXI Horizontal Unit

With Genaxy Scientific GENECO range of gel electrophoresis units, you can perform analytical & preparative separations of nucleic acids simply & efficiently. Each system has a wide variety of comb types & comes with a typically waterproof gel casting system. To fulfill your needs, we've created electrophoresis cells that are remarkably advanced, versatile & user-friendly.



**GENECO-H-MAXI**

## Technical Specification

Model	GENECO-H-MAXI
Principal Material	Acrylic
Dimensions (L X B X H)	35 X 22 X 8 cm
No. of Samples	26
Gel Size	20 x 20 cm
Trays (L x W)	20 X 20 cm (1)
Combs	26 Well x 1.5mm (1) Acrylic 3 Well Preparative x 3 mm (1) Acrylic 3 Well Preparative x 6 mm (1) Acrylic
Gel Casting Tray	1 Universal
Platinum Electrode	Red & Black, 1 each
Connecting Cord	Red & Black, 1 each
Lid	1 No.

## Horizontal System Includes

Basic Unit (1 No.)

Platinum Electrode Assembly Removable (2 Nos.)

Connecting Cord (1 Nos.)

26 well combs 1.5 mm (2 Nos.)

Gel Running Tray 200 X 200 mm (1 No.)

Lid (1 No.)

3 Well Preparative Comb 3.0 mm (1 No.)

3 Well Preparative Comb 6.0 mm (1 No.)

### Accessories

S No.	Catalog No.	Description
1	GENECO-H-MAXI-20	Gel Running tray 20 x 20 cm
2	GENECO-H-MAXI-20UV	UV transparent Gel Running tray 20 x 20cm
3	GENECO-H-LT-25	Levelling table for gel casting 25 x 25 cm
4	GENECO-H-MAXI-EL	Spare electrode 2 Nos.

### Combs

S No.	Code	Description
1	GENECO-H-MAXI-1.5-26	Acrylic Comb MINI 26 Well 1.5mm
2	GENECO-H-MAXI-6.0-3	Acrylic Comb MINI 3 Well 6 mm
3	GENECO-H-MAXI-3.0-3	Acrylic Comb MINI 3 Well 3 mm
4	GENECO-H-MAXI-CS	Comb Stand

## Genaxy Scientific Pvt. Ltd

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com