

WIDE BORE TIPS



They're a perfect fit for popular pipettors like Eppendorf, Gilson, and Capp. Wide bore tips, also known as large orifice tips, are the go-to choice in labs for tackling thick or dense liquids, like proteins, gels, and genomic DNA. They're especially handy for pipetting high-viscosity samples, working with cell cultures, and handling protein and DNA tasks. These tips are essential for the routine processing of samples, making them a practical choice for a variety of laboratory applications.

- Available Sizes: 200uL, 1000uL.
- Packaging Options: Ready in racks or bags.

The wide orifice helps keep delicate samples safe by lowering the force, great for dealing with thick stuff.

Cat. Number	Description	Qty./Pack
GEN-UT-205-WB-C	200 μ L, Edgeless Large-bore, Natural, Bulk	1000 Tips
GEN-UT-205-WB-R-S	200 μ L, Edgeless Large-bore, Natural, Racked, Sterile	96x10 Racks
GEN-UT-1005-WB-C	1000 μ L, Edgeless Large-bore, Natural	1000 Tips
GEN-UT-1005-WB-R-S	1000 μ L, Edgeless Large-bore, Natural, Racked, Sterile	96x10 Racks

Genaxy Scientific Pvt. Ltd.