



“ just keep spinning ”



CENTRIFUGE TUBES

ELITE PERFORMANCE

— AT —

EXCEPTIONAL VALUE

GENAXY Brand Centrifuge Tubes are designed to meet the needs Of Researchers. This clear polypropylene Centrifuge tubes comes with the advanced High Density Polyethylene caps packed in high quality ziplock bags. These tubes are available in 7ml, 10ml, 15 ml, 25ml, 30 ml and 50 ml volume with amber (15ml, 50ml Tubes Only) and clear color options. Genaxy also provide these tubes in sterile, non-sterile, racked and individually wrapped formats.


Unique Features

- » Premium grade polypropylene resin conforming to US FDA 21 CFR
- » Manufactured Under ISO Class 8 Clean Room Facility
- » High spin rates with maximum RCF ratings of up to 20,000 x g for 50ml Tubes and 17,000 x g for 15ml Tubes
- » Temperature ranges from -80°C to 121°C
- » Easy Grip Caps have an integral seal ring design for secure closures that eliminate seepage when used as recommended
- » Clear, Black, accurate printed graduations with unique Eco-friendly Hot Stamp Technology
- » Cap can be closed with quick 2-3 turns allows single hand operation for closing and opening of the tube.
- » Large white writing area for sample labelling.
- » Ultra-smooth conical bottoms, which allow for 100% visibility and no sharp "gates" to cause cuts and scratches.
- » Sterile tubes are radiation sterilized which ensures the non-cytotoxicity.
- » Certified DNase/RNase, Heavy-metal, Human DNA, Phthalate and pyrogen free
- » Gamma irradiated to SAL 10-6



Genaxy Centrifuge Tubes

ECO-FRIENDLY HOT STAMP PRINTING



Flat Cap For Labeling


Easy Grip Cap with integral seal ring for easy gripping and one hand operation



Writing Area For Sample Identification



Extremely Accurate Printed Graduations With Eco-Friendly Hot Stamp Printing



Ultra Smooth Bottom For 100% Visibility- No Sharp 'Gates' To Cause Cuts And Scratches



FOR MORE PRODUCT INFORMATION
SCAN THE QR CODE





Cat. Number	Description	Qty. / Pack
GEN-CT-15ML	15ml Centrifuge Tube with Cap, Conical, Clear, Non-Sterile	25 X 20
GEN-CT-15ML-S	15ml Centrifuge Tube with Cap, Conical, Clear, Sterile	25 X 20
GEN-CT-15ML-R	15ml Centrifuge Tube with Cap, Conical, Racked, Clear, Non-Sterile	25 X 20
GEN-CT-15ML-R-S	15ml Centrifuge Tube with Cap, Conical, Racked, Clear, Sterile	25 X 20
GEN-CT-25ML-S	25ml Centrifuge Tube with Cap, Conical, Clear, Sterile	25 X 8
GEN-CT-25ML-R-S	25ml Centrifuge Tube with Cap, Conical, Racked, Clear, Non-Sterile	25 X 8
GEN-CT-50ML	50ml Centrifuge Tube with Cap, Conical, Clear, Non-Sterile	25 X 20
GEN-CT-50ML-S	50ml Centrifuge Tube with Cap, Conical, Clear, Pre- Sterilized	25 X 20
GEN-CT-15ML-X	15ml Centrifuge Tube with Cap, Conical, Amber, Non-Sterile	25 X 20
GEN-CT-15ML-X-S	15ml Centrifuge Tube with Cap Conical, Amber, Sterile	25 X 20
GEN-CT-15ML-X-R	15ml Centrifuge Tube with Cap, Conical, Racked, Clear, Non-Sterile	25 X 20
GEN-CT-15ML-X-R-S	15ml Centrifuge Tube with Cap, Conical, Racked, Clear, Sterile	25 X 20
GEN-CT-50ML-X	50ml Centrifuge Tube with Cap, Conical, Amber, Non-Sterile	25 X 20
GEN-CT-50ML-X-S	50ml Centrifuge Tube with Cap, Conical, Amber, Sterile	25 X 20
GEN-CT-50ML-R	50ml Centrifuge Tube with Cap, Conical, Racked, Clear, Non-Sterile	15 X 20
GEN-CT-50ML-R-S	50ml Centrifuge Tube with Cap, Conical, Racked, Clear, Sterile	15 X 20
GEN-CT-50ML-SS	50ml Centrifuge Tube with Cap, Self Standing, Clear, Non-Sterile	25 X 20
GEN-CT-50ML-SS-S	50ml Centrifuge Tube with Cap, Self Standing, Clear, Sterile	25 X 20
GEN-CT-15ML-S-IW	15ml Centrifuge Tube with Cap, Conical, Clear, Individually Wrapped, Sterile	250
GEN-CT-50ML-S-IW	50ml Centrifuge Tube with Cap, Conical, Clear, Individually Wrapped, Sterile	250
GEN-CT-30ML-SS-S	30ml Centrifuge Tube with Cap, Self Standing, Clear, Sterile	500
GEN-CT-7ML-U	7ml Centrifuge Tube with Plug Cap, Round Bottom, Clear	2000
GEN-CT-10ML-U	10ml Centrifuge Tube with Plug Cap, Round Bottom, Clear	2000
GEN-CT-50ML-U	50ml Centrifuge Tube with Plug Cap, Round Bottom, Clear	500

Genaxy Scientific Pvt. Ltd.