

GENAXY Scientific Pvt. Ltd. is a leading manufacturer of laboratory glassware and plasticware. We offer a wide range of products including microcentrifuge tubes, PCR tubes, and pipette tips. Our products are designed for high precision and low retention, making them ideal for a variety of laboratory applications. For more information, please visit our website at www.genaxy.com.

UNIVERSAL FIT GENTIPS

LOW RETENTION

HIGH PRECISION


GENTip universal tip comes with a bevelled or fine point orifice. Bevelled orifice which gives shorter contact area while aspirating and dispensing, results in No Last Drop Issue and No False Results. While the standard tips, without bevelled orifice, comes with larger surface area, which reduces the control, results in last drop issues and no accuracy. GenTips are available in bulk packs, racked form and individually packed. Bulk tips are bagged, and heat sealed in puncture resistant laminated nylon-base dispenser packs with super wide gusset at bottom to stand neatly upright. Stack rack options are also available in 10ul, 200ul and 1000ul.

Unique Features

- » Made of 99.9% virgin Polypropylene.
- » Manufactured under ISO Class 8 clean room facility
- » Bevelled orifice for total sample recovery
- » Unique unigrip design with 3 internal rings creating a vacuum allowed onto the pipette shaft eliminating sideways dispensing
- » Use of high-quality diamond polished infection mould for smooth internal surface and reduced liquid hang-up
- » Zero Fluid Retention (ZFR) technology offers maximum fluid recovery.
- » Accurate graduation marks provide visual confirmation that the correct liquid volume is being aspirated or dispensed
- » Touch centrum design allows tips to fit all the single and multichannel pipettors
- » Tips are DNase/RNase, Human DNA, Heavy metal, Pyrogen, Phthalate and PCR inhibitor free.
- » USFDA 21 CFR Compliant



Genaxy Universal Fit GenTips



Beveled Orifices For Total Sample Recovery

Blade Prevents Droplets & Increases Precision



Molded In Exclusive 99.9% Virgin Polypropylene Resin Which Results In Exceptional Sample Delivery & Improved Assay Performance

High-quality diamond polished injection mould for smooth internal surface and reduced liquid hang-up



Better Sealability With Soft But Touch Centrum Design



FOR MORE PRODUCT INFORMATION
SCAN THE QR CODE



Universal Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-UT-10-C	0.5-10µl Universal Fit Gentip, Natural, Bulk, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-10-L-C	0.5-10µl Universal Fit Gentip, Natural, Bulk, Low Retention, Non-Sterile, 1000 tips/Pk	1000 tips
GEN-UT-10-C-R	0.5-10µl Universal Fit Gentip Natural, Racked, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-10-C-R-S	0.5-10µl Universal Fit Gentip, Natural, Racked, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-10-C-STK	0.5-10µl Universal Fit Gentip, Natural, Reload Stacks, Non-Sterile, 96 tips/Stack, 960 tips/Pk	960 tips
GEN-UT-10-L-C-STK	0.5-10µl Universal Fit Gentip, Natural, Reload Stacks, Low Retention, Non-Sterile, 96 tips/Stack, 960 tips/Pk,	960 tips
GEN-UT-10-C-L-R	0.5-10µl Universal Fit Gentip, Natural, Racked, Low Retention, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-10-C-L-R-S	0.5-10µl Universal Fit Gentip, Natural, Racked, Low Retention, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-10XL-C	0.5-20µl Universal Fit Gentip, Natural, Bulk, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-10XL-L-C	0.5-20µl Universal Fit Gentip, Natural, Bulk, Low Retention, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-10XL-C-R	0.5-20µl Universal Fit Gentip, Natural, Racked, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-10XL-C-R-S	0.5-20µl Universal Fit Gentip, Natural, Racked, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-10XL-C-L-R	0.5-20µl Universal Fit Gentip, Natural, Racked, Low Retention, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-10XL-C-L-R-S	0.5-20µl Universal Fit Gentip, Natural, Racked, Low Retention, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-200-C	1-200µl Universal Fit Gentip, Natural, Bulk, Bevelled, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-200-Y	1-200µl Universal Fit Gentip, Yellow, Bulk, Bevelled, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-200-C-NB	1-200µl Universal Fit Gentip, Natural, Bulk, Non-Bevelled, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-200-Y-NB	1-200µl Universal Fit Gentip, Yellow, Bulk, Non-Bevelled, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-200-L-C	1-200µl Universal Fit Gentip, Natural, Bulk, Bevelled, Low Retention, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-200-C-R	1-200µl Universal Fit Gentip, Natural, Racked, Graduated, Bevelled, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks

Universal Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-UT-200-Y-R	1-200µl Universal Fit Gentip, Yellow, Racked, Bevelled, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-200-C-R-S	1-200µl Universal Fit Gentip, Natural, Racked, Graduated, Bevelled, Sterile., 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-200-Y-R-S	1-200µl Universal Fit Gentip, Yellow, Racked, Bevelled, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-200-C-L-R	1-200µl Universal Fit Gentip, Natural, Racked, Graduated, Bevelled, Low Retention, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-200-C-L-R-S	1-200µl Universal Fit Gentip, Natural, Racked, Graduated, Bevelled, Low Retention, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-200-C-STK	1-200µl Universal Fit Gentip, Natural, Reload Stacks, Bevelled, Non-Sterile, 96 tips/Stack, 960 tips/Pk	960 tips
GEN-UT-200-Y-STK	1-200µl Universal Fit Gentip, Yellow, Reload Stacks, Bevelled, Non-Sterile, 96 tips/Stack, 960 tips/Pk	960 tips
GEN-UT-200-LC-STK	1-200µl Universal Fit Gentip, Natural, Reload Stacks, Bevelled, Low Retention, Non-Sterile, 96 tips/Stack, 960 tips/Pk	960 tips
GEN-UTG-200-C	1-200µl Universal Fit Gentip, Natural, Bulk, Graduated, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UTG-200-Y	1-200µl Universal Fit Gentip, Yellow, Bulk, Graduated, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UTG-200-L-C	1-200µl Universal Fit Gentip, Natural, Bulk, Graduated, Low Retention, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UTG-200-C-R	1-200µl Universal Fit Gentip, Natural, Racked, Graduated, Bevelled, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-200-Y-R	1-200µl Universal Fit Gentip, Yellow, Racked, Bevelled, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-200-C-R-S	1-200µl Universal Fit Gentip, Natural, Racked, Graduated, Bevelled, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-200-Y-R-S	1-200µl Universal Fit Gentip, Yellow, Racked, Bevelled, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-200-C-L-R	1-200µl Universal Fit Gentip, Natural, Racked, Graduated, Bevelled, Low Retention, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-200-C-L-R-S	1-200µl Universal Fit Gentip, Natural, Racked, Graduated, Bevelled, Low Retention, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-200XL-C	1-200µl (89mm) Universal Fit Gentip, Natural, Bulk, Graduated, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UTG-200XL-C-R	1-200µl (89mm) Universal Fit Gentip, Natural, Racked, Graduated, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-200XL-C-R-S	1-200µl (89mm) Universal Fit Gentip, Natural, Racked, Graduated, Sterile, 96 tips/rack, 10 racks/pk	10 racks

Universal Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-UT-300-C	1-300µl Universal Fit Gentip, Natural, Bulk, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-300-L-C	1-300µl Universal Fit Gentip, Natural, Bulk, Low Retention, Non-Sterile, 1000 tips/pk	1000 tips
GEN-UT-300-C-R	1-300µl Universal Fit Gentip, Natural, Racked Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-300-C-R-S	1-300µl Universal Fit Gentip, Natural Racked, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-300-C-L-R	1-300µl Universal Fit Gentip, Natural, Racked, Bevelled, Low Retention, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-300-C-L-R-S	1-300µl Universal Fit Gentip, Natural, Racked, Bevelled, Low Retention, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UT-1000-C	100-1000µl Universal Fit Gentip, Natural, Bulk, Bevelled, Non-Sterile, 500 tips/pk	500 tips
GEN-UT-1000-B	100-1000µl Universal Fit Gentip, Blue, Bulk, Bevelled, Non-Sterile, 500 tips/pk	500 tips
GEN-UT-1000-L-C	100-1000µl Universal Fit Gentip, Natural, Bulk, Bevelled, Low Retention, Non-Sterile, 500 tips/pk	500 tips
GEN-UT-1000-C-R	100-1000µl Universal Fit Gentip, Natural, Racked, Bevelled, Non-Sterile, 100 tips/rack, 10 racks/pk	10 racks
GEN-UT-1000-B-R	100-1000µl Universal Fit Gentip, Blue, Racked, Bevelled, Non-Sterile, 100 tips/rack, 10 racks/pk	10 racks
GEN-UT-1000-C-R-S	100-1000µl Universal Fit Gentip, Natural, Racked, Bevelled, Sterile, 100 tips/rack, 10 racks/pk	10 racks
GEN-UT-1000-B-R-S	100-1000µl Universal Fit Gentip, Blue, Racked, Bevelled, Sterile, 100 tips/rack, 10 racks/pk	10 racks
GEN-UT-1000-C-L-R	100-1000µl Universal Fit Gentip, Natural, Racked, Bevelled, Low Retention, Non-Sterile, 100 tips/rack, 10 racks/pk	10 racks
GEN-UT-1000-C-L-R-S	100-1000µl Universal Fit Gentip, Natural, Racked, Bevelled, Low Retention, Sterile, 100 tips/rack, 10 racks/pk	10 racks
GEN-UTG-1000-C	1-1000µl Universal Fit Gentip, Natural, Graduated, Bulk, Non-Sterile, 500 tips/pk	500 tips
GEN-UTG-1000-B	1-1000µl Universal Fit Gentip, Blue, Graduated, Bulk, Non-Sterile, 500 tips/pk	500 tips
GEN-UTG-1000-L-C	1-1000µl Universal Fit Gentip, Blue, Graduated, Low Retention, Bulk, Non-Sterile, 500 tips/pk	500 tips
GEN-UTG-1000-C-STK	1-1000µl Universal Fit Gentip, Natural, Reload Stacks, Graduated, Non-Sterile, 96 tips/stack, 960 tips/pk	960 tips
GEN-UTG-1000-LC-STK	1-1000µl Universal Fit Gentip, Natural, Reload Stacks, Graduated, Low Retention, Non-Sterile, 96 tips/stack, 960 tips/pk	960 tips

Universal Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-UTG-1000-C-R	100-1000µl Universal Fit Gentip, Natural, Racked, Graduated, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-1000-B-R	100-1000µl Universal Fit Gentip, Blue, Racked, Graduated, Non-Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-1000-C-R-S	100-1000µl Universal Fit Gentip, Natural, Racked, Graduated, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-1000-B-R-S	100-1000µl Universal Fit Gentip, Blue, Racked, Graduated, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-1000-C-L-R	100-1000µl Universal Fit Gentip, Natural, Racked, Graduated, Low Retention, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-1000-C-L-R-S	100-1000µl Universal Fit Gentip, Natural, Racked, Graduated, Low Retention, Sterile, 96 tips/rack, 10 racks/pk	10 racks
GEN-UTG-1250-C	100-1250µl Universal Fit Gentip, Natural, Bulk, Graduated, Non-Sterile (104mm), 500 tips/pk	500 tips
GEN-UTG-1250-B	100-1250µl Universal Fit Gentip, Blue, Bulk, Graduated, Non-Sterile (104mm), 500 tips/pk	500 tips
GEN-UTG-1250-L-C	100-1250µl Universal Fit Gentip, Natural, Bulk, Graduated, Low Retention, Non-Sterile (104mm), 500 tips/pk	500 tips
GEN-UTG-1250-C-R	100-1250µl Universal Fit Gentip, Natural, Racked, Graduated, Non-Sterile.(104mm), 96 tips/rack, 8 racks/pk	8 racks
GEN-UTG-1250-B-R	100-1250µl Universal Fit Gentip, Blue, Racked, Graduated, Non-Sterile.(104mm), 96 tips/rack, 8 racks/pk	8 racks
GEN-UTG-1250-B-R-S	100-1250µl Universal Fit Gentip, Blue, Racked, Graduated, Sterilized. (104mm), 96 tips/rack, 8 racks/pk	8 racks
GEN-UTG-1250-C-R-S	100-1250µl Universal Fit Gentip, Natural, Racked, Graduated, Sterilized.(104mm), 96 tips/rack, 8 racks/pk	8 racks
GEN-UTG-1250-C-L-R	100-1000µl Universal Fit Gentip, Natural, Racked, Bevelled, Retention, Sterile, 96 tips/rack, 8 racks/pk	8 racks
GEN-UTG-1250-C-L-R-S	100-1250µl Universal Fit Gentip, Natural, Racked, Graduated, Low Retention, Sterile.(104mm), , 96 tips/rack, 8 racks/pk	8 racks
GEN-UT-1250XL-C	100-1200µl Universal Fit Gentip, Natural, Bulk, Non-Sterile (120mm), 500 tips/pk	500 tips
GEN-UT-1250XL-B	100-1200µl Universal Fit Gentip, Blue, Bulk, Non-Sterile (120mm), 500 tips/pk	500 tips
GEN-UT-1250XL-L-C	100-1200µl Universal Fit Gentip, Natural, Bulk, Non-Sterile (120mm), Low Retention, 500 tips/pk	500 tips
GEN-UT-1250XL-C-R	100-1200µl Universal Fit Gentip, Natural, Racked, Non-Sterile (120mm), 96 tips/rack, 6 rack/pk	6 racks
GEN-UT-1250XL-B-R	100-1200µl Universal Fit Gentip, Blue, Racked, Non-Sterile.(120mm), 96 tips/rack, 6 rack/pk	6 racks

Cat. Number	Description	Qty. / Pack
GEN-UT-1250XL-C-R-S	100-1200µl Universal Fit Gentip, Natural, Racked, Sterile.(120mm), 96 tips/rack, 6 rack/pk	6 racks
GEN-UT-1250XL-B-R-S	100-1200µl Universal Fit Gentip, Blue, Racked, Sterile.(120mm), 96 tips/rack, 6 rack/pk	6 racks
GEN-UT-1250XL-C-L-R	100-1200µl Universal Fit Gentip, Natural, Racked, (120mm), Low Retention, 96 tips/rack, 6 rack/pk	6 racks
GEN-UT-1250XL-C-L-R-S	100-1200µl Universal Fit Gentip, Natural, Racked, (120mm), Low Retention, Sterile, 96 tips/rack, 6 rack/pk	6 racks

Empty Racks for Pipette Tips

GEN-UT-10-ER	Empty Racks for GenTip 10µl	10 racks
GEN-UT-200-ER	Empty Racks for GenTip 200µl	10 racks
GEN-UT-300-ER	Empty Racks for GenTip 300µl	10 racks
GEN-UT-1000-ER	Empty Racks for GenTip 1000µl	10 racks
GEN-UTG-1000-ER	Empty Racks for GenTip 1200µl	10 racks
GEN-UTG-1250-ER	Empty Racks for GenTip 1250µl	8 racks
GEN-UTG-1250XL-ER	Empty Racks for GenTip 1250µl XL	6 racks

Individually Wrapped Pipette Tips

GEN-C-300-IW	10µl Universal Fit Gentip, Clear, Individual Peel Off, Plastic Wrap, Sterile, 500 tips/pk	500 tips
GEN-Y-200-IW	200µl Universal Fit Gentip, Yellow, Individual Peel Off, Plastic Wrap, Sterile, 500 tips/pk	500 tips
GEN-C-200-IW	200µl Universal Fit Gentip, Clear, Individual Peel Off Plastic Wrap, Sterile, 500 tips/pk	500 tips
GEN-C-200-IW-PP	200µl Universal Fit Gentip, Clear, Individual Peel Off, Paper Plastic Wrap, Sterile, 2000 tips/pk	5 x 400 tips
GEN-B-1000-IW	1000µl Universal Fit Gentip, Blue, Individual Peel Off, Plastic Wrap, Sterile, 500 tips/pk	500 tips
GEN-C-1000-IW	1000µl Universal Fit Gentip, Clear, Individual Peel Off, Plastic Wrap, Sterile, 500 tips/pk	500 tips
GEN-C-1000-IW-PP	1000µl Universal Fit Gentip, Clear, Individual Peel Off, Paper Plastic Wrap, Sterile, 1000 tips/pk	1000 tips
GEN-C-5000-IW	5000µl Universal Fit Gentip, Clear, Individual Peel Off, Plastic Wrap, Sterile, 100 tips/pk	100 tips
GEN-C-5000-IW-PP	5000µl Universal Fit Gentip, Clear, Individual Peel Off, Paper Plastic Wrap, Sterile, 100 tips/pk	100 tips

Cat. Number	Description	Qty. / Pack
GEN-UT-5000-C	1-5ml Universal Fit Gentip, Natural, Bulk, Gilson Style, Non-Sterile, 100 tips/pk	100 tips
GEN-UT-5000-R-S	1-5ml Universal Fit Gentip, Natural, Racked, Gilson Style, Sterile, 24 tips/rack, 8 racks/pk	8 racks
GEN-GPL-5WB-C	5mL Pipette Tip, Wide Bore, 130mm length, bagged, For Gilson PL, 100 tips/pk	100 tips
GEN-GPL-5WB-R-S	5mL Pipette Tip, Wide Bore, 130mm length, Racked, For Gilson PL, 24 tips/rack, 8 racks/pk	8 racks
GEN-ET-5WB-C	Long 5mL Pipette Tip, Wide Bore, 175mm length, bagged, For Eppendorf/Sartorius, 100 tips/pk	100 tips
GEN-ET-5WB-R-S	Long 5mL Pipette Tip, Wide Bore, 175mm length, racked, For Eppendorf/Sartorius, 24 tips/rack, 8 racks/pk	8 racks
GEN-TT-5SB-C	5mL Pipette Tip, Small Bore, 150mm length, bagged, For Brand/Thermo, 100 tips/pk	100 tips
GEN-TT-5SB-R-S	5mL Pipette Tip, Small Bore, 150mm length, racked, For Brand/Thermo, 40 tips/rack, 8 racks/pk	8 racks
GEN-ETS-5WB-C	Short 5mL Pipette Tip, Wide Bore, 130mm length, bagged, For Eppendorf/Sartorius/GilsonP, 100 tips/pk	100 tips
GEN-ETS-5WB-R-S	Short 5mL Pipette Tip, Wide Bore, 130mm length, racked, For Eppendorf/Sartorius/GilsonP, 24 tips/rack, 8 racks/pk	8 racks
GEN-UT-10ML-C	1-10ml Universal Fit Gentip, Natural, Bulk, Gilson Style, Non-Sterile, 200 tips/pk	200 tips
GEN-UT-10ML-R-S	1-10ml Universal Fit Gentip, Natural, Racked, Gilson Style, Non-Sterile, 24 tips/rack, 8 racks/pk	8 racks
GEN-TT-10WB-C	10mL Pipette Tip, Wide Bore, 161mm length, bagged, For Thermo, 100 tips/pk	100 tips
GEN-TT-10WB-R-S	10mL Pipette Tip, Wide Bore, 161mm length, racked, For Thermo, 24 tips/rack, 8 racks/pk	8 racks
GEN-ET-10SB-C	10mL Pipette Tip, Small Bore, 165mm length, bagged, For Eppendorf/Sartorius/Brand/GilsonP,PN,PL,PG/RaininUNV, 100 tips/pk	100 tips
GEN-ET-10SB-R-S	10mL Pipette Tip, Small Bore, 165mm length, racked, For Eppendorf/Sartorius/Brand/GilsonP,PN,PL,PG/RaininUNV, 24 tips/rack, 8 racks/pk,	8 racks

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com