

GENAXY PAPER ELECTROPHORESIS

With Genaxy Scientific GENECO range of gel electrophoresis units, you can perform analytical & preparative separations of nucleic acids simply & efficiently. It has a wide variety of comb types & comes with a typically waterproof gel casting system. To fulfill your needs, we've created electrophoresis cells that are remarkably advanced, versatile & user-friendly.

Basic unit consists of anode and cathode buffer reservoirs capable of accommodating suitable electrode assembly. In case of vertical model paper supporting vertical stands with cross rods are provided in between Reservoirs, whereas a removable horizontal supporting stand has been Provided in horizontal model in both the cases transparent acrylic cover are provided for controlling buffer evaporation and for electrical safety.

Catalog No.	Description	Quantity
GENECO-PE-VMINI	Paper Electrophoresis System (Vertical Mini for 3 Strips)	01 Unit
GENECO-PE-VLARGE	Paper Electrophoresis System (Vertical Large for 6 Strips)	01 Unit
GENECO-PE-HLARGE	Paper Electrophoresis System (Horizontal Large for 6 strips)	01 Unit

**PAPER ELECTROPHORESIS SYSTEM****Genaxy Scientific Pvt. Ltd**

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com