



## ROD/DISC GEL ELECTROPHORESIS SYSTEMS

	GENECO-ROD-RDual	GENECO-ROD-CS	GENECO-ROD-CL
<b>Principal Material</b>	Acrylic	Acrylic	Acrylic
<b>Dimensions</b>	13 X 13 X 7cm	10 X 14 cm (id X h)	15 X 24 cm (id X h)
<b>Buffer Volume</b>	250 + 250 mL	250 mL	500 + 500 ml
<b>Glass tubes</b>	14 cm long X 5mm i.d (12 nos)	7.5cm long X 5mm i.d (8 nos)	14cm long X 5mm i.d (12 nos)
	7.5cm long X 5mm i.d (12 nos)		7.5 cm long X 5mm i.d (12 nos)
<b>Rubber Grommets</b>	12 Nos.	8 Nos.	12 Nos.
<b>Tube Stand Number</b>	1 No	1 No	1 No
<b>Tube Stand Type</b>	Adjustable	Fixed	Adjustable
<b>Long Needles</b>	2 Nos.	2 Nos.	2 Nos.
<b>No. Of Platinum Electrode</b>	Red and Black (1 Each)	Red and Black (1 Each)	Red and Black (1 Each)
<b>Connecting Cord</b>	Red and Black (1 Each)	Red and Black (1 Each)	Red and Black (1 Each)
<b>Lid</b>	1 No	3 No	1 No
<b>Gel Tubes</b>		8 Tubes	12 Tubes



**GENAXY SCIENTIFIC PVT. LTD.**

155-156, C-Block, Plot No. 32-57, Community Centre, Behind Janak Cinema,  
Janak Puri, New Delhi-110058, Phone:011-4761999 (100 lines), Fax:011-25534163, Email:sales@genaxy.com