



GENAXY VERTICAL SYSTEMS

MAXI Vertical Unit

With Genaxy Scientific GENECO range of gel electrophoresis units, you can perform analytical & preparative separations of nucleic acids simply & efficiently. Each system has a wide variety of comb types & comes with a typically waterproof gel casting system. To fulfill your needs, we've created electrophoresis cells that are remarkably advanced, versatile & user-friendly.

GENECO-V Vertical Systems



GENECO offers a complete range of Vertical Gel Units that combines the ease of use with high resolution separations. REGULAR as well as DUAL Models are available to run two gels simultaneously.

Features:

- Compatible with all type of Precast gels of major brands
- Rapid gel casting & loading
- Low buffer volumes required
- Main unit is of Acrylic & Combs & spacers are made up of Teflon

Technical Specification

Model	GENECO-V-MAXI	GENECO-V-MAXI DUAL
Principal Material	Acrylic	Acrylic
Dimensions (L X B X H)	29 X 20 X 21 cm	29 X 22 X 23 cm
No. of Samples	20	40
Gel Size	16 x 20 cm	16 x 20 cm (2 Gel)
Buffer Volume	300 + 600 ml	300 + 700 ml
Combs	20 Well – Teflon 0.5 mm, 1 mm, 1.5 mm (1)	20 Well – Teflon 0.5 mm, 1 mm (2)
Teflon Spacers	0.5 mm, 1 mm, 1.5 mm (6)	0.5 mm, 1 mm (4)
Platinum Electrode	Red & Black, 1 each	Red & Black, 1 each
Connecting Cord	Red & Black, 1 each	Red & Black, 1 each
Lid	1 No.	1 No.
Glass Plates	Notch/Rectangular (2) Set	Notch/Rectangular (2) Set
Gasket	Fixed	Fixed
Clamp/Screws	2 Sets	4 Sets
Water Circulation	No	Yes

Accessories

S No.	Catalog No.	Description
1	GENECO-V-MAXI-GP/NP	Glass Plates (Notched or Regular)
2	GENECO-V-MAXI-Caster	Gel Caster for 2 Gels
3	GENECO-V-MAXI-GDT	Gel Developing Tray with gel holder
4	GENECO-V-MAXI-EL	Spare Electrode 2 Nos.
5	GENECO-V-MAXI-GPS	Glass Plates Stand for Maxi System

Combs & Spacers

S No.	Code	Description
1	GENECO-V-MAXI- x-20	Teflon Combs MAXI 20 Well
2	GENECO-V-MAXI- x	Teflon Spacers (6 Nos/Pack)

Genaxy Scientific Pvt. Ltd

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com