


SINGLE CHANNEL PIPETTE GENsleek


Variable Volume



CALIBRATION AND ADJUSTMENT

Take out the Service tool and place the Service tool as shown in the figure, Turn the service tool clockwise to increase, or counterclockwise to decrease the volume.

Calibrate pipettes 1-2 times per year
(depending on frequency of use)



Product Features

- Single channel and Multi-Channel are available.
- Lightweight and Ergonomic design provides excellent operating, very suitable for gripping. Light spring, press very labour-saving, help to relieve fatigue of hand joints, easy suction;
- Digital display, adjusted delivery volume clearly indicated in the large digital display on the handle;
- Fully made of high strength material. Anti-strengthening corrosion, alcohol disinfection; Enhanced anti UV linearity, UV sterilization;
- Different colors are used to identify pressing handles, making it easier to select pipettes of required specifications;
- Well-compatible with most of the pipette tips on the market;
- The pipette is equipped with GenSleek pipette tips to ensure good air tightness and accurate pipetting volume.

Cat. No.	Description	Volume/Increment	Qty./Pk	Price/Pk
GENSleek-2	Digital Micro Pipette	0.1-2.5 μ l / 0.002 μ l	each	6,630.00
GENSleek-10	Digital Micro Pipette	0.5-10 μ l / 0.1 μ l	each	6,630.00
GENSleek-20	Digital Micro Pipette	2-20 μ l / 0.5 μ l	each	6,630.00
GENSleek-50	Digital Micro Pipette	5-50 μ l / 0.5 μ l	each	6,630.00
GENSleek-100	Digital Micro Pipette	10-100 μ l / 1 μ l	each	6,630.00
GENSleek-200	Digital Micro Pipette	20-200 μ l / 1 μ l	each	6,630.00
GENSleek-1000	Digital Micro Pipette	100-1000 μ l / 5 μ l	each	6,630.00
GENSleek-5000*	Digital Micro Pipette	500-5000 μ l / 50 μ l	each	6,630.00
GENSleek-10000*	Digital Micro Pipette	1000-10000 μ l / 0.1ml	each	6,630.00
GENSleek-T-5ML	Tips for GenSleek-5000, Clear	1-5ml	100 tips	1,990.00
GENSleek-T-10ML	Tips for GenSleek-10000, Clear	1-10ml	100 tips	2,330.00
GEN-SAS-01	Carousel Stand for 6 Single Channel Pipettes	-	each	3,735.00

*Order Specialized tips from Genaxy