



CRYOGENIC VIALS

*Where Samples Find
Eternal Preservation*

Genaxy Cryogenic Vials are intended to store samples in freezers or similar low-temperature equipment. During storage, build up of stress on the product can happen, the genaxy design as well as the high grade material ensures safe handling during and after the unfreezing process.

External Thread, Cryo Tubes

Features

- » Without/O-rings with the geometry of the external thread, Genaxy can abstain from using a silicone gasket,
- » The external threading results in a smooth internal surface that eliminates undue loss of specimen material during removal.
- » Supplied in convenient, puncture and tear resistant re-closable bags with zipper and laser perforation
- » Self-standing, Star foot
- » Graduated with large writing surface.

Internal Thread, Cryo Tubes with 'O' ring

Features

- » Made from special polypropylene (PP),
- » These 2.0ml (working volume 1.8ml) cryovials can withstand temperatures ranging from 121°C to -196 °C.
- » The self-standing tubes have internal threading with maximum graduation of 1.80ml



Cryogenic Vial

Cat. Number	Description	Qty. / Pack
GEN-CV-1.2ml-E	1.2ml Cryogenic Vial, External Thread, Graduated, Sterile.	50 x 10
GEN-CV-1.8ml-E	1.8ml Cryogenic Vial, External Thread, Graduated, Sterile.	50 x 10
GEN-CV-1.8ml-EX	1.8ml Cryogenic Vial, External Thread, Graduated, Amber, Sterile.	50 x 10
GEN-CV-5.0ml-E	5.0ml Cryogenic Vial, External Thread, Graduated, Sterile.	25X10
GEN-CVE-COLOR	Cap Inserts for External Thread Vials, Color: Red, Violet, Green, Yellow, Blue	100
GEN-CV-1.2ml-I	1.2ml Cryogenic Vial, Internal Thread, Graduated, Sterile.	50x10
GEN-CV-1.8ml-I	1.8ml Cryogenic Vial, Internal Thread, Graduated, Sterile.	50x10
GEN-CV-1.8ml-IR	1.8ml Cryogenic Vial, Round Bottom, Internal Thread, Graduated, Sterile.	50x10
GEN-CV-1.8ml-IX	1.8ml Cryogenic Vial, Internal Thread, Graduated, Amber, Sterile.	50x10
GEN-CV-5.0ml-I	5.0ml Cryogenic Vial, Internal Thread, Graduated, Sterile.	25X10
GEN-CVI-COLOR	Cap Inserts for Internal Thread Vials, Color: Red, Violet, Green, Yellow, Blue	100

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
 Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com