



MINI TUBE SYSTEM


INDIVIDUALLY UNIQUE

TOGETHER COMPLETE


The MicroTube Rack System is the industry standard for deep-well microtiter plates with removable wells. These racks can be used for DNA elution in automatic system in high throughput isolation, individual elution, DNA storage and transportation etc. Good product for Seed Industries, Biotechnology Labs, Genomic Labs etc.

Unique Features

- » Mini Tube System is made up of 99.9% virgin polypropylene
- » Manufactured under ISO Class 8 clean room facility
- » SBS compatible, 96 tube slots
- » Available as single tubes or in strips of 8 tubes
- » Compatible with robotic workstations (Like Beckman, Tecan, etc.), multi-channel pipettors and auto-sampling instruments.
- » An alpha-numeric grid, non-reversing rack and lid ensure positive sample identification.
- » The lid, with individual condensation rings above each tube, prevents cross-contamination.
- » The inert polypropylene tubes, rack and lid can be autoclaved and centrifuged at moderate speeds
- » Sterilized as well as non-sterilized tubes, strip tubes, tube racks and strip caps are available
- » Can be stacked and used for long term storage of biological samples
- » Certified DNase, RNase, Human DNA and Pyrogen free
- » US FDA 21 CFR compliant



Genaxy Mini Tube System



An alpha-numeric grid, non-reversing rack and lid ensure positive sample identification.

8 X12 hole pattern on traditional 9mm centers



Strip of 8 tubes



Individual Tubes



FOR MORE PRODUCT INFORMATION SCAN THE QR CODE



Mini Tube System

Cat. Number	Description	Qty. / Pack
GEN-MT-65-C	0.65ml Mini Tube, Single Tubes, Natural, Bulk, Non-Sterile	4800
GEN-MT-65-C-R	0.65ml Mini Tube, Single Tubes, Natural, Racked (96 Tube-Slot), Non-Sterile	4800
GEN-MT-65-C-R-S	0.65ml Mini Tube, Single Tubes, Natural, Racked (96 Tube-Slot), Sterile	4800
GEN-MT-6508-C	0.65ml Mini Tube, 8-Strips Tubes, Natural, Bulk, Non-Sterile	600 strips
GEN-MT-01-C	1.1ml Mini Tube, Single Tubes, Natural, Bulk, Non-Sterile	4800
GEN-MT-01-C-R	1.1ml Mini Tube, Single Tubes, Natural, Racked (96 Tube-Slot), Non-Sterile	4800
GEN-MT-01-C-R-S	1.1ml Mini Tube, Single Tubes, Natural, Racked (96 Tube-Slot), Sterile	4800
GEN-MT-0108-C	1.1ml Mini Tube, 8-Strips Tubes, Natural, Bulk, Non-Sterile	600 strips
GEN-MT-0108-C-R	1.1ml Mini Tube, 8-Strips Tubes, Natural, Racked (96 Tube-Slot), Non-Sterile	4800
GEN-MT-0108-C-R-S	1.1ml Mini Tube, 8-Strips Tubes, Natural, Racked (96 Tube-Slot), Sterile	4800
GEN-MTC-8-NVC	8-Strips Mini Tube Caps, PE, Non-Vented, Natural, Non-Sterile	600 strips
GEN-MTC-8-NVC-S	8-Strips Mini Tube Caps, PE, Non-Vented, Natural, Sterile	600 strips

Genaxy Scientific Pvt. Ltd.