

# PS600 Power Supply (600 Volt)

The versatile PS600 is designed for most electrophoresis techniques

## Convenient and Simple

The small footprint, large handle, and simple operation makes the PS600 easy to set up and use in your lab. With the large, easy to read LED, simply set the constant parameter, set the timer (if desired), and push start.



## Specifications

<b>Power requirement</b>	Frequency, 50–60 Hz Line voltage, 90–260 VAC
<b>User interface</b>	3 digit, red LED membrane keypad
<b>Output</b>	5–600 V, 1 V increment 3–750 mA, 1 mA increment 150 W maximum 4 pairs output jacks in parallel
<b>Output connector</b>	4 mm plugs
<b>Timer</b>	1–999 min, 1 min increment
<b>Environmental operating conditions</b>	Indoor use: 0–40 °C Relative humidity: 0–80% noncondensing Altitude: < 2000 m Pollution degree: 2 Installation category: II
<b>Dimensions (h × w × d)</b>	33.5 × 12.5 × 28 cm 13.2 × 4.9 × 11.0 in
<b>Weight</b>	2.4 kg (5.3 lb)
<b>Safety protection circuits</b>	Short circuit protection Open circuit protection
<b>Product certifications</b>	EN61010-1:2010, EN61326-1:2006 Class B CE WEEE RoHS