

## Specimen or Sample Transport Box

Ideal for transporting biological and clinical samples, as well as instruments and products that you want to keep clean and dry under testing environmental conditions

The DuraPorter® transparent design is OSHA 29 CFR 1910.1031 and UN3373 compliant

- Polycarbonate construction and available in three colors for user/owner identification
- Accommodates 13 mm, 16 mm and 20 mm tubes with our Half and Full Size OneRack®
- Use the separators provided to create your own compartment sizes



Item No.	Color		L x W x H in	L x W x H cm
120052	Clear with Blue Handle	■		
120077	Red with Red Handle	■	15 x 7.7 x 6.3	38 x 19.6 x 16
120078	Yellow with Yellow Handle	■		





Load



Lock



Transport



**Store**  
can be placed in a freezer or refrigerator to keep samples chilled

**Lid Opens a Full 180°**  
to expedite loading and emptying





**Extra Security Feature**  
Accommodates lockout tags (not included)



Included separators can also be used to create inner compartments.  
Or, use without racks or separators as one large storage space for larger samples or supplies

**A Perfect Fit**  
**OneRack® Tube Rack**  
Accommodates 13 mm, 16 mm and 20 mm tubes with our Half and Full Size OneRack®



Corporate Office :

**GENAXY SCIENTIFIC PVT. LTD.**

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058

Ph. : +91-11-4761 9999 (100Lines) Web. : [www.genaxy.com](http://www.genaxy.com) E-mail : [technical@genaxy.com](mailto:technical@genaxy.com)