

Biohazard Disposal Bags

For personal or general lab use, sizes to suit all needs

- Polypropylene
- Includes an indicator patch that darkens when autoclaved



Epoxy Coated
HDPE Wire Holder

Biohazard Disposal Bags

Item No.	Color	Thickness mil	Capacity	Packaging	L x W in	L x W cm
HS10320			2 L	500	8 x 12	20.3 x 30.5
HS10321	Red ■	1.6	9 L	200	14 x 19	35.6 x 48.3
HS10322			19 L	200	19 x 23	48.3 x 58.4
HS10323			40 L	200	25 x 33	63.5 x 83.8

HDPE Coated Steel Wire Bag Holder

Item No.	L x W x H in	L x W x H cm
HS10044H	3.6 x 5 x 8.4	9.1 x 12.8 x 21.3

*Holds Biohazard Disposal Bag HS10320



Autoclavable Disposal Bags

Clear bags imprinted with blue biohazard labeling designed for safe disposal

- Constructed of 2 mil thick polypropylene, allowing for repeated autoclaving
- Front write-on panel allowing for repeated autoclaving
- Strong side seams for securely holding items
- 4 inner packs of 50 bags



Case of 200

Item No.	Color	Capacity	L x W in	L x W cm
HS1002A		14 L	12.2 x 26	31 x 66
HS1002B	Clear □	26 L	16.3 x 23.6	41.5 x 60
HS1002C		69 L	24 x 31.9	62 x 81

Corporate Office :

GENAXY SCIENTIFIC PVT. LTD.

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058

Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com