

### GENAXY Pipette Controller

- Can fill 25ml pipet under 4 seconds
- Safety valve and hydrophobic filters provide double protection against fluid penetration
- Two different speed modes-High & Low along with Gravity drain
- The pump speed can also be fine-tuned by varying finger pressure on operating knobs for better control of speed
- Cadmium free environment friendly NiMH batteries
- Batteries can be changed very easily by the user
- The intelligent charger prevents over charging/heating of batteries
- Low battery indicator
- The specially designed Desktop Stand enables the Pipet Controller to be charged while resting on it.
- The Pipet Controller can hold the serological pipette while resting on Desktop stand.



Catalog No	Description	Volume
GENPET-PM	Pipette Controller with Soft grip push buttons, Double safety valve and Hydrophobic PTFE Membrane Filter 0.45 $\mu$ m & 0.20 $\mu$ m	1-100 ml

### Genaxy Scientific Pvt. Ltd

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com