

PIPETTE DISPENSER

It is a lightweight (approx. 174g) and slim pipette controller that is easy to grip.

Discharge can be switched between natural discharge and forced discharge by pressing the discharge button, enabling seamless discharge operation.

It uses a lithium-ion battery that enables long-term operation, and can operate for about 8 hours with intermittent suction and discharge per full charge.

A wall mount stand is included.

Specifications

Product Code	GENPET-epf
Product Name	Pipette Mate Slim
Suction/Discharge	The suction speed can be changed to 6 levels, and the discharge speed is constant (natural discharge and forced discharge can be switched by pressing or lowering the button)
filter	φ25mm, pore size 0.45μm, hydrophobic
battery	Lithium-ion rechargeable battery (non-removable, please contact your dealer if you wish to replace it)
USB-AC adapter for charging	Input: AC100~240V, 50/60Hz, Output: DC5V 1.0A (Please be sure to use the included USB-AC adapter for charging)
Standard Accessories	Filter (installed inside and 1 spare included), wall mount stand

