


# Clinical 4-Way Tube Rack


The only tube rack specifically designed to hold round, flat & conical bottom blood collection and culture tubes

## Product Highlights

- Accommodates all brands of 13 and 16 mm round, flat and conical bottom blood collection & culture tubes
- Six way connection, link different or the same well sizes, end to end or side by side, in any orientation
- Customizable racks can be linked together, or used as a standalone
- Assorted color pack available for color coding options
- Racks can be used in water baths and placed in freezers



Holds  
Multiple Tube Sizes  
without multiple  
racks!



**Assorted Color Pack**  
available for color coding options

### Blood Tube Brands (include & not limited to):


- Becton Dickinson
- Greiner
- Covidien
- Thermo Scientific™
- Applied Biosystems
- Sarstedt

### Recommended Tube Dimensions:

- 13 x 75 mm
- 13 x 100 mm
- 16 x 100 mm
- 16 x 125 mm
- 17 x 120 mm


### Tube Volume Guide:

- 3.5 mL
- 7.0 mL
- 10 mL
- 15 mL



### Customizable

Racks can slide together to create an infinite number of combinations




## Clinical 4-Way Tube Rack

- Accommodates all brands of 13 and 16 mm round, flat and conical bottom blood collection & culture tubes
- Only 4-way rack available with a six way connection, link different or the same well sizes, end to end or side by side, in any orientation
- Customize the rack to fit the number of wells required, infinite units can be linked together or use one standalone unit
- Uniformly spaced wells help eliminate pipetting errors which enables more accurate results
- Molded alphanumeric identification for quick location and retrieval of tubes
- Robust, made from polypropylene that is chemical resistant and easy to clean
- Racks can be used in water baths and placed in freezers
- Assorted color pack available for color coding options
- Patent Pending
- **Autoclavable**



**Molded Alphanumeric Identification**  
for quick location and retrieval of tubes



**Tube Bottoms Nest**  
in well pockets keeping them from going thru the rack and off the bench top



### Versatile


Accommodates tubes with round, flat and conical bottoms

↑  
**Holds 15 mL Conical Tubes!**

### Individual Rack Specifications

Item No.	Color	Item No.	Color	Side	Wells	Well Dia. mm	L x W x H in	cm	UOM
120684	Blue	120687	Natural	A	12	13	9.1 x 2.8 x 3.8	23.2 x 7.0 x 9.5	5 ea/pk, 6 pk/cs
120685	Purple	120683	Assorted*	B	18	13			
120686	Green			C	12	17			
				D	18	17			

\*Assorted colors (Blue, Purple, Green, Red, Natural)



Corporate Office

**GENAXY SCIENTIFIC PVT. LTD.**

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058

Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com