

Micro Weigh Boats

Ideal for microbalance weighing and/or distribution of very small amounts of sampling material


Easy to Handle

- Large tabs make boats easy to use, position and hold, especially with forceps

Superb Sample Recovery

- Inert material enables the boat and contents to be added to solutions to aid 100% sample recovery, where protocols allow
- Anti-static properties prevent sample retention from static electricity
- Small surface area reduces liquid sample evaporation

Actual size



Volume
500 µl




Easy to Handle Tabs

Made of Polystyrene

- Food Grade Compliant

Holds approx 500 µl



Ideal for very small volume pipette calibration and checks and micro sampling

Micro Weigh Boats

Pack of 250

Item No.	Color	Volume µl	L x W x H in	L x W x H mm
120836	White	500	.98 x 0.4 x 0.3	25 x 8.9 x 7.5

GENAXYSCIENTIFIC PVT. TD.

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058
Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com