


X-tourn® Tourniquets

Hi-viz orange alert for safer practice

High-quality disposable tourniquets eliminate roping or curling when used

- Made from supple and smooth polyisoprene
- Latex free
- Less abrasive than other brands
- No rubber proteins, resins or fatty acids that can produce hypersensitive reactions
- Available in bulk or on a perforated roll



Uniform Thickness

prevents shredding and breaking

| Item No. | Color | Packaging | L x W in | L x W cm |
|----------|--|-----------|-------------|-------------|
| 120108 | Green ■ | 100/Bulk | 18 x 1 | 45.7 x 2.5 |
| 120110 | | 250/Bulk | | |
| 120111 | Orange ■ | 25 Rolls | 18 x 1 | 45.7 x 2.5 |