


Genaxy Gel Loading Tips



Gel Loading tips are designed for easy loading of samples on gels. The tips fit a wide variety of single and multichannel pipettors. Tips are available in flat and round formats with outer diameter dimensions for specific gel thickness. All tips are certified DNase, RNase and Pyrogen free.

Unique Features

- » Microflex Gel Loading Tips are designed for easy loading of samples on gels
- » Available in flat and round formats with outer dimensions for specific gel thicknesses.
- » Tips comes in three different outer dimensions : 0.17mm, 0.4mm and 0.57mm
- » Fits into wide variety of single and multichannel pipettors.
- » Certified DNase-/RNase-free and nonpyrogenic

Cat. Number	Description	Qty. / Pack
GEN-TGL-120R	10µl MicroFlex Gel Loading Tips, Natural, Racked, Flat Tip 0.17mm, Non-Sterile	4 X 200 Tips
GEN-TGL-120R-S	10µl MicroFlex Gel Loading Tips, Natural, Racked, Flat Tip 0.17mm, Sterile	4 X 200 Tips
GEN-TGL-320R	200µl MicroFlex Gel Loading Tips, Natural, Racked, Flat Tip 0.4mm, Non-Sterile	4 X 200 Tips
GEN-TGL-320R-S	200µl MicroFlex Gel Loading Tips, Natural, Racked, Flat Tip 0.4mm, Sterile	4 X 200 Tips
GEN-TGL-520R	200µl MicroFlex Gel Loading Tips, Natural, Racked, Round Tip 0.57mm, Non-Sterile	10 X 96 Tips
GEN-TGL-520R-S	200µl MicroFlex Gel Loading Tips, Natural, Racked Round Tip 0.57mm, Sterile	10 X 96 Tips

Genaxy Scientific Pvt. Ltd.