



TRANSPORT TUBES

MOVE SAFELY


— *WITH* —

US

Genaxy durable and reliable Transport Tube and Cap provide a 'leak-proof' and reliable seal for your critical sample storage. Both Tube and Cap are manufactured of virgin polypropylene giving our product the ability to be autoclaved.

Unique Features

- » Manufactured with High grade 99.9% Virgin polypropylene under ISO class 8 clean room facility
- » Leak-proof and reliable seal for critical sample storage
- » Graduation marks at 1 ml intervals
- » Large, frosted writing area for sample labelling
- » Sterile or Non-sterile forms are available
- » Sterile Tubes are either ETO sterile or Gamma sterile (both options are given)
- » 5ml and 12ml volumes are available
- » Certified DNase/RNase, Human DNA and Phthalate free
- » US FDA 21 CFR Compliant



Genaxy Transport Tubes



Leak-proof and reliable seal ensures safe transportation of samples



Accurate graduation for visual confirmation



Frosted area for writing



Skirted bottom for better stability on benchtop



FOR MORE PRODUCT INFORMATION
SCAN THE QR CODE



Transport Tubes

| Cat. Number | Description | Qty. / Pack |
|----------------|--|-------------|
| GEN-ST-5ml | 5ml Transport Tubes, Clear, Non-Sterile | 500 |
| GEN-ST-5ml-S | 5ml Transport Tubes, Clear, Sterile (ETO Sterile) | 500 |
| GEN-ST-5ml-GS | 5ml Transport Tubes, Clear, Sterile (GAMMA Sterile) | 500 |
| GEN-ST-12ml | 12ml Transport Tubes, Clear, Non-Sterile | 500 |
| GEN-ST-12ml-S | 12ml Transport Tubes, Clear, Sterile | 500 |
| GEN-ST-12ml-GS | 12ml Transport Tubes, Clear, Sterile (GAMMA Sterile) | 500 |

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com