

GENAXY ROUND BOTTOM TEST TUBE



17x100mm, 14 mL polystyrene round bottom test tubes have a 1,400 RCF rating. The tubes come with polyethylene dual-position snap caps, which offer both vented and fully closed positions. RNase/DNase-free.

Widely used and referenced in laboratory protocols
Dual position snap-cap prevents sample loss with ease of handling

Sterile, Non-pyrogenicity tested to less than 0.1 EU/ ml

Cat. Number	Description	Volume	Size (ØxH) mm	Type/Style	Qty. / Pack
GEN-14100-PS	Test Tubes, Round Bottom	14ml	17 X 100	PS	500
GEN-14100-PS-IW	Test Tubes, Round Bottom, Individually Wrapped	14ml	17 X 100	PS	500

Genaxy Scientific Pvt. Ltd.