



DEEP WELL PLATES

*AN IDEAL SOLUTION
FOR*

HIGH THROUGHPUT RESEARCH


Deepwell Plates are an ideal product for short and long term sample storage, tissue culture work, enzyme assays, and extraction of nucleic acids and proteins and screening studies. Genaxy deep well plates are conform to SBS standards for use with multichannel pipettes and automatic systems.

Unique Features


- » Made from USPVI Polypropylene material
- » Conforms to SBS standards
- » Compatible with multichannel pipettes and automatic equipment.
- » High chemical stability.
- » Temperature ranges from -80 °C to 120 °C.
- » Centrifugal strength goes up to 3,000 xg .
- » 96 -Deep-Well Plate are available with Round Conical/U-bottom
- » Optimized bottom design to minimize residual volume.
- » Clear alphanumeric codes for easy identification.
- » Stackable
- » DNase / RNase, Human DNA and pyrogen free



Genaxy Deep Well Plates




U-Bottom



Conical V-Bottom



Elution Plate



96-Well Comb



8-Well Comb



FOR MORE PRODUCT
INFORMATION
SCAN THE QR CODE



Deep Well Plates

Cat. Number	Description	Qty. / Pack
GEN-96-0.5-RC	200ul, 96-Deep-Well Plate square well Round Conical bottom (Elution plate), Suitable with Thermo System (KingFisher Flex)	50
GEN-96-1.6-SC	1.6ml, 96-Deep-Well Plate square well Round Conical bottom (Elution plate), Suitable with Genetix System	50
GEN-96-2.2-SC	2.2ml 96-Deep-Well Plate square well Round Conical bottom, Suitable with Thermo System (KingFisher Flex) and Genetix System	50
GEN-96-2.2-SU	2.2ml 96-Deep-Well Plate square well Round U-bottom, Suitable with Zybio System	50
GEN-YMC-96	96 Magnetic tip comb, Suitable with Thermo System (KingFisher Flex)	100/Pack
GEN-YMC-8	8 Magnetic tip comb, Suitable with Zybio System	2 Pcs/Bag, 240 Bags/Cs

Genaxy Scientific Pvt. Ltd.