

## PARATAPE WAXFILM



WAXFILM is a polymer lipid film made of paraffin composite materials. It is a thermoplastic with extensibility, waterproofing, odorlessness, translucency, and cohesiveness, and has versatility under various conditions and aromatic solvents. The drum of Waxfilm is made of resin material, so it will not produce debris and affect the experimental data. Waxfilm has better stretchability and is less prone to breakage. Waxfilm has better outer Layer Paper and do not break easily when wet. The drum of Waxfilm is made of resin material, so it will not produce debris.

Cat. Number	Description	Dimensions	Qty. / Pack
GEN-PSF-001	ParaTape Sealing Film, 4 IN. x 125 FT	10 cm x 38.1 mtr.	1 Roll
GEN-PSF-002	ParaTape Sealing Film, 2 IN.x 165 FT	5 cm x 50 mtr.	1 Roll
GEN-PSF-003	ParaTape Sealing Film, 1.2 IN.x 50 FT	3 cm x 15 mtr.	1 Roll
GEN-PSF-004	ParaTape Sealing Film, 0.4 IN.x 50 FT	1 cm x 15 mtr.	1 Roll

## Genaxy Scientific Pvt. Ltd.