

SEROLOGICAL PIPETTES



High accuracy disposable serological pipettes in polystyrene. One piece construction, Colour coded plug, Fine graduations increase accuracy and ease reading. Permanent and relief printing, Accuracy >98% at full volume, pipettes suited for use with any pipetting device, Pipettes are gamma irradiated, Sterile pipettes are non pyrogenic, non cytotoxic, non haemolytic, comes with Certificates of conformity

- Made of high clarity, 100% virgin polystyrene.
- Sterilized by E-beam, SAL=10⁻⁶.
- Non-Pyrogenic, DNase/Rnase free.
- Edged color-coded band for identification.
- Negative graduations for extra volume.
- Polyolefin fiber filter to reduce contamination.

Cat. Number	Description	Volume	Qty. / Pack
GEN-SER-1ML-S	Individually wrapped Serological Pipette, Sterile	1ml	500
GEN-SER-2ML-S	Individually wrapped Serological Pipette, Sterile	2ml	400
GEN-SER-5ML-S	Individually wrapped Serological Pipette, Sterile	5ml	200
GEN-SER-10ML-S	Individually wrapped Serological Pipette, Sterile	10ml	200
GEN-SER-25ML-S	Individually wrapped Serological Pipette, Sterile	25ml	200
GEN-SER-50ML-S	Individually wrapped Serological Pipette, Sterile	50ml	100

Genaxy Scientific Pvt. Ltd.