

SCREW CAP TUBES

*BEING VERSATILE
IS A STYLE*


GENAXY screw cap micro-tubes serve a wide range of laboratory needs including collection, processing, centrifugation, long term storage and packaging/shipping of restriction enzymes and other reagents. Conical & self standing both options are offered with 7 different color alternatives (White, Pink, Blue, Green, Orange, Red and Yellow) for caps.

Unique Features

- » Manufactured with High grade 99.9% Virgin polypropylene.
- » Manufactured Under ISO class 8 cleanroom facility
- » Suitable for High-Speed Centrifugation up to 20,000xg
- » Can be used at temperature ranging from -195 deg C (LN2 vapour phase) to 121 deg C
- » The threaded "O" ring design, containing Silicon rubber, ensures a uniform closure
- » Outer Walls of the Tube are ribbed which aids in single handed operation of the same
- » Screw Cap Tubes Are Now Also Available With Looped Caps
- » These unique caps have a circular loop which makes closing and opening of the tube easier and prevents dropping caps on the laboratory floor
- » Caps have extra long lips to help protect samples
- » Certified DNase/RNase and Human DNA and phthalate free
- » USFDA 21CFR Compliant



Genaxy Screw Cap Tubes



These unique caps have a circular loop which makes closing and opening of the tube easier and prevents dropping caps on the laboratory floor

The threaded "O" ring design to ensures uniform closure

Outer walls of the tube are ribbed for single handed operation

Ultra smooth tube bottoms for 100% visibility - No Sharp Gates to cause cuts and scratches



FOR MORE PRODUCT INFORMATION SCAN THE QR CODE



Screw Cap Tubes

Cat. Number	Description	Qty. / Pack
GEN-ST-050-C	0.5ml Screw Cap Tubes with Caps, Conical, Natural, Non-Sterile	500 X 2
GEN-ST-050-C-S	0.5ml Screw Cap Tubes with Caps, Conical, Natural, Sterile	10 X 100
GEN-ST-050-A	0.5ml Screw Cap Tubes with Caps, Conical, Assorted*, Non-Sterile	500 X 2
GEN-ST-050-A-S	0.5ml Screw Cap Tubes with Caps, Conical, Assorted*, Sterile	10 X 100
GEN-ST-050-COLOR	0.5ml Screw Cap Tubes with Caps, Conical, W/P/R/G/Y/B, Non-Sterile	500 X 2
GEN-ST-050-X	0.5ml Screw Cap Amber Tubes with Caps, Conical, Brown, Non-Sterile	500 X 2
GEN-ST-050-X-S	0.5ml Screw Cap Amber Tubes with Caps, Conical, Brown, Sterile	10 X 100
GEN-ST-050-SS-C	0.5ml Screw Cap Tubes with Caps, Self Standing, Natural, Non-Sterile	500 X 2
GEN-ST-050-SS-C-S	0.5ml Screw Cap Tubes with Caps, Self Standing, Natural, Sterile	10 X 100
GEN-ST-050-SS-A	0.5ml Screw Cap Tubes with Caps, Self Standing, Assorted, Non-Sterile	500 X 2
GEN-ST-050-SS-A-S	0.5ml Screw Cap Tubes with Caps, Self Standing, Assorted, Sterile	10 X 100
GEN-ST-050-SS-COLOR	0.5ml Screw Cap Tubes with Caps, Self Standing, W/P/R/G/Y/B, Non-Sterile	500 X 2
GEN-ST-050-SS-X	0.5ml Screw Cap Amber Tubes with Caps, Self Standing, Brown, Non-Sterile	500 X 2
GEN-ST-050-SS-X-S	0.5ml Screw Cap Amber Tubes with Caps, Self Standing, Brown, Sterile	10 X 100
GEN-ST-045-SS	0.5ml Elongated Type Screw Cap Vial, Self Standing, Non-Sterile	500 X 2
GEN-ST-045-SS-S	0.5ml Elongated Type Screw Cap Vial, Self Standing, Sterile	10 X 100
GEN-ST-150-C	1.5ml Screw Cap Tubes with Caps, Conical, Natural, Non-Sterile	500 X 2
GEN-ST-150-C-S	1.5ml Screw Cap Tubes with Caps, Conical, Natural, Sterile	10 X 100
GEN-ST-150-A	1.5ml Screw Cap Tubes with Caps, Conical, Assorted*, Non-Sterile	500 X 2
GEN-ST-150-A-S	1.5ml Screw Cap Tubes with Caps, Conical, Assorted*, Sterile	10 X 100
GEN-ST-150-COLOR	1.5ml Screw Cap Tubes with Caps, Conical, W/P/R/G/Y/B, Non-Sterile	500 X 2
GEN-ST-150-X	1.5ml Screw Cap Amber Tubes with Caps, Conical, Brown, Non-Sterile	500 X 2
GEN-ST-150-X-S	1.5ml Screw Cap Amber Tubes with Caps, Conical, Brown, Sterile	10 X 100
GEN-ST-150-SS-C	1.5ml Screw Cap Tubes with Caps, Self Standing, Natural, Non-Sterile	500 X 2
GEN-ST-150-SS-C-S	1.5ml Screw Cap Tubes with Caps Self Standing, Natural, Sterile	10 X 100
GEN-ST-150-SS-A	1.5ml Screw Cap Tubes with Caps Self Standing, Assorted, Non-Sterile	500 X 2
GEN-ST-150-SS-A-S	1.5ml Screw Cap Tubes with Caps Self Standing, Assorted, Sterile	10 X 100
GEN-ST-150-SS-COLOR	1.5ml Screw Cap Tubes with Caps, Self Standing, W/P/R/G/Y/B, Non-Sterile	500 X 2
GEN-ST-150-SS-X	1.5ml Screw Cap Amber Tubes with Caps, Self Standing, Brown, Non-Sterile	500 X 2
GEN-ST-150-SS-X-S	1.5ml Screw Cap Amber Tubes with Caps, Self Standing, Brown, Sterile	10 X 100
GEN-ST-200-C	2.0ml Screw Cap Tubes with Caps, Conical, Natural, Non-Sterile	500 X 2
GEN-ST-200-C-S	2.0ml Screw Cap Tubes with Caps, Conical, Natural, Sterile	10 X 100
GEN-ST-200-A	2.0ml Screw Cap Tubes with Caps, Conical, Assorted*, Non-Sterile	500 X 2
GEN-ST-200-A-S	2.0ml Screw Cap Tubes with Caps, Conical, Assorted*, Sterile	500 X 2
GEN-ST-200-COLOR	2.0ml Screw Cap Tubes with Caps, Conical, W/P/R/G/Y/B, Non-Sterile	500 X 2
GEN-ST-200-X	2.0ml Screw Cap Amber Tubes with Caps, Conical, Brown, Non-Sterile	500 X 2
GEN-ST-200-X-S	2.0ml Screw Cap Amber Tubes with Caps, Conical, Brown, Sterile	500 X 2

Assorted* _ White, Pink, Blue, Green, Orange, Red, Yellow

Screw Cap Tubes

Cat. Number	Description	Qty. / Pack
GEN-ST-200-SS-C	2.0ml Screw Cap Tubes with Caps, Self Standing, Natural, Non-Sterile	500 X 2
GEN-ST-200-SS-C-S	2.0ml Screw Cap Tubes with Caps, Self Standing, Natural, Sterile	10 X 100
GEN-ST-200-SS-A	2.0ml Screw Cap Tubes with Caps, Self Standing, Assorted, Non-Sterile	500 X 2
GEN-ST-200-SS-A-S	2.0ml Screw Cap Tubes with Caps, Self Standing, Assorted, Sterile	10 X 100
GEN-ST-200-SS-COLOR	2.0ml Screw Cap Tubes with Caps, Self Standing, W/P/R/G/Y/B,	500 X 2
	Non-Sterile	
GEN-ST-200-SS-X	2.0ml Screw Cap Amber Tubes with Caps, Self Standing, Brown,	500 X 2
	Non-Sterile	
GEN-ST-200-SS-X-S	2.0ml Screw Cap Amber Tubes with Caps, Self Standing, Brown, Sterile	10 X 100

Screw Cap Micro Tubes with Looped Caps

Cat. Number	Description	Qty. / Pack
GEN-ST-050-C-LP	0.5ml Screw Cap Tubes with Looped Caps (with O-ring) , Conical, Natural, Non-Sterile, Unassembled	500 X 2
GEN-ST-050-C-S-LP	0.5ml Screw Cap Tubes with Looped Caps (with O-ring) , Conical, Natural, Sterile, Assembled	10 X 100
GEN-ST-050-SS-C-LP	0.5ml Screw Cap Tubes with Looped Caps (with O-ring) , Self Standing, Natural, Non-Sterile, Unassembled	500 X 2
GEN-ST-050-SS-C-S-LP	0.5ml Screw Cap Tubes with Looped Caps (with O-ring) , Self Standing, Natural, Sterile, Assembled	10 X 100
GEN-ST-045-SS-LP	0.5ml Elongated Type Screw Looped Cap Vial with Caps (with O-ring) , Self Standing, Non-Sterile., Unassembled	500 X 2
GEN-ST-045-SS-S-LP	0.5ml Elongated Type Screw Looped Cap Vial with caps (with O-ring) , Self Standing, Sterile, Assembled	10 X 100
GEN-ST-150-C-LP	1.5ml Screw Cap Tubes with Looped Caps (with O-ring) , Conical, Natural, Non-Sterile, Unassembled	500 X 2
GEN-ST-150-C-S-LP	1.5ml Screw Cap Tubes with Looped Caps (with O-ring) , Conical, Natural, Sterile, Assembled	10 X 100
GEN-ST-150-SS-C-LP	1.5ml Screw Cap Tubes with Looped Caps (with O-ring) , Self Standing, Natural, Non-Sterile, Unassembled	500 X 2
GEN-ST-150-SS-C-S-LP	1.5ml Screw Cap Tubes with Looped Caps (with O-ring) Self Standing, Natural, Sterile, Assembled	10 X 100
GEN-ST-200-C-LP	2.0ml Screw Cap Tubes with Looped Caps (with O-ring) , Conical, Natural, Non-Sterile, Unassembled	500 X 2
GEN-ST-200-C-S-LP	2.0ml Screw Cap Tubes with Looped Caps (with O-ring) , Conical, Natural, Sterile, Assembled	10 X 100
GEN-ST-200-SS-C-LP	2.0ml Screw Cap Tubes with Looped Caps (with O-ring) , Self Standing, Natural, Non-Sterile, Unassembled	500 X 2
GEN-ST-200-SS-C-S-LP	2.0ml Screw Cap Tubes with Looped Caps (with O-ring) , Self Standing, Natural, Sterile, Assembled	10 X 100
GEN-SCO-A-LP	Looped Screw Caps (with O-ring) Color* (W/P/R/G/Y/B)	500 X 2

Genaxy Scientific Pvt. Ltd.