

REAGENT RESERVOIR



- "Trough within a trough" maximizes reagent recovery
- Convenient pour-off spouts on all four corners
- A wide base provides stability and avoids spills
- Graduated marks on inside wall
- Made from high impact polystyrene
- Superior design provides stronger sidewalls
- Sterile product is sterilized by gamma irradiation,
- RNASE and DNASE free
- 25ml divider allows for pipetting with eight channels on one side and four channels or a single channel on the other; the divider minimizes reagent loss and is convenient for adding controls

Cat. Number	Description	Volume	Color	Sterility	Qty. / Pack
GEN-RR-25	Reagent reservoir	25ml	White	Non-Sterile	100
GEN-RR-25D	Reagent reservoir with divider	25ml	White	Non-Sterile	100
GEN-RR-50	Reagent reservoir	50ml	White	Sterile	100
GEN-RR-100	Reagent reservoir	100ml	White	Non-Sterile	100

Genaxy Scientific Pvt. Ltd.