



Premium Cryo Tubes



External Thread, Cryo Tubes W/O-rings

These Cryogenic vials are manufactured from polypropylene to withstand temperatures to -196°C (-321°F). All vials can be color-coded with cap inserts, which are sold separately. Vials have graduation marks and a large marking spot for easy labeling of contents. Supplied Sterile. DNase and RNase-free, and nonpyrogenic.

Warning : Do not use cryo vials in the liquid phase of liquid nitrogen.

Catalog No.	Description	Capacity (ml)	Qty/Pack
GEN-CVE-1.2ml	Cryogenic Vial, 1.2ml, Sterilized, Graduated	1.2ml	50 X 10 Tubes
GEN-CVE-2.0ml	Cryogenic Vial, 2.0ml, Sterilized, Graduated	2.0ml	50 X 10 Tubes
GEN-CVE-5.0ml	Cryogenic Vial, 5.0ml, Sterilized, Graduated	5.0ml	25 X 10 Tubes
GEN-CVE-1.2ml-X	Cryogenic Vial, 1.2ml, Sterilized, Graduated, Amber Color	1.2ml	50 X 10 Tubes
GEN-CVE-2.0ml-X	Cryogenic Vial, 2.0ml, Sterilized, Graduated, Amber Color	2.0ml	50 X 10 Tubes
GEN-CVE-5.0ml-X	Cryogenic Vial, 5.0ml, Sterilized, Graduated, Amber Color	5.0ml	25 X 10 Tubes

Internal Thread, Cryo Tubes W/O-rings

These Cryogenic vials are manufactured from polypropylene to withstand temperatures to -196°C (-321°F). All vials can be color-coded with cap inserts, which are sold separately. Vials have graduation marks and a large marking spot for easy labeling of contents. Supplied Sterile. DNase- and RNase-free, and nonpyrogenic.

Warning : Do not use cryo vials in the liquid phase of liquid nitrogen.

Catalog No.	Description	Capacity (ml)	Qty/Pack
GEN-CVI-1.2ml	Cryogenic Vial, 1.2ml, Sterilized, Graduated	1.2ml	50 X 10 Tubes
GEN-CVI-2.0ml	Cryogenic Vial, 2.0ml, Sterilized, Graduated	2.0ml	50 X 10 Tubes
GEN-CVI-5.0ml	Cryogenic Vial, 5.0ml, Sterilized, Graduated	5.0ml	25 X 10 Tubes
GEN-CVI-1.2ml-X	Cryogenic Vial, 1.2ml, Sterilized, Graduated, Amber Color	1.2ml	50 X 10 Tubes
GEN-CVI-2.0ml-X	Cryogenic Vial, 2.0ml, Sterilized, Graduated, Amber Color	2.0ml	50 X 10 Tubes
GEN-CVI-5.0ml-X	Cryogenic Vial, 5.0ml, Sterilized, Graduated, Amber Color	5.0ml	25 X 10 Tubes

Pickup Device for Cryo Vials

Catalog No.	Description	Qty/Pack
GEN-PD-IN	Pickup Device for Internal Thread Cryo Vials	10/Pk
GEN-PD-EX	Pickup Device for External Thread Cryo Vials	10/Pk

Cap Inserts for Internally Thread Vials

Catalog No.	Description	Qty/Pack
GEN-BC-CI-Y	Cap Inserts for Internal Thread Vials, Yellow	50/Pk
GEN-BC-CI-L	Cap Inserts for Internal Thread Vials, Lime Green	50/Pk
GEN-BC-CI-P	Cap Inserts for Internal Thread Vials, Pink	50/Pk
GEN-BC-CI-PU	Cap Inserts for Internal Thread Vials, Purple	50/Pk