

## Packed with features for comfortable and accurate pipetting

The powerful pump offers rapid, yet sensitive aspiration and dispensing, has adjustable valves and is microprocessor controlled for constant reliability and accuracy

- Long-life li-ion battery, with quick charge—up to 20 hours use
- Toggle operation by pressing the large mode button
- An adapter can be fitted to the nose cone to stabilize small volume and thinner pipettes
- New USB connector ensures the appropriate charger is used with the unit
- Low-battery indicator
- Can be recharged when in use



Rechargeable li-ion battery (replaceable)



Nose cone assembly only

Item No.	Volume Range mL	L x W x H in	L x W x H cm
HS3000	1–100	9 x 8.5 x 5	22.9 x 21.6 x 12.7

### Replacement parts

HS5000RFFIL	RF3000® Replacement Filters, Pack of 5		
120669	RF3000® Universal Power Adaptor 5VDC Micro USB Plug		
HS300502	RF3000® Li-Ion Rechargeable Battery, One year warranty on battery included		
120487	RF3000® Nose Cones Assorted* ■ ■ ■ ■ ■ Set of 5		
HS5000RFRFA	RF3000® Replacement Filter Adapter		
HS5000RFRP	RF3000® Replacement Pipette Silicone Adapter		

\*Assorted colors (Black, Blue, Grey, Red and Yellow)

### Specifications

<b>Power</b>	Micro-USB connector Rechargeable li-ion battery (replaceable)
<b>Filter</b>	PTFE/25mm dia/0.45 µm filter
<b>Pipettes</b>	1–100 mL, glass and plastic
<b>Pipette Nozzle</b>	Polypropylene

### Package Includes

RF3000® Pipette Controller
Charging Stand and Wall-Mountable/Magnetic Holder
Nozzle Adapter for use with Small Pipettes
Li-Ion Rechargeable Battery
Five Nose Cones: Black, Yellow, Red, Grey and Blue
5VDC Micro USB Plug
6 ft (183 cm) Cable

1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



**Universal Grip**  
for right- or left-handed users

**Soft-Coated Triggers**  
with direction indicators

**Ergonomically Balanced**  
in the hand with a loaded pipette

**UV resistant housing**

**Benchtop charging base**



### Easy Identification

5 different colored nose cones help differentiate department, task or user; as well as a labeling area on the side of the unit

**Adapter**  
for small pipettes

**Storage Stand**  
doubles as a charging station



**Magnetic & Screw Wall-Mount**



**3 Modes**  
Bright LEDs indicate 3 operation modes: high, low, & gravity

Corporate Office  
**GENAXY SCIENTIFIC PVT. LTD.**

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058  
Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com