


smartPicks® Inoculation Picks

The 'smart' way to isolate individual colonies

Available in three lengths and two needle sizes to suit a wide range of application needs

- Naturally biodegradable, made from 100% U.S.A. grown white birch wood
- Manufactured in an FDA-registered factory and complies with FDA standards for food contact
- Available in semi or needle point styles for delicate operations



Mini Sized SmartPicks

fit neatly into a 1.5 mL microfuge tube allowing the lid to be closed for full sample transfer/recovery

Drop it in & Close it up!

Item No.	Size	Point	Tip	Suggested Tube mL	Packaging	Length in	Length cm
120010	Mini	Needle Point		Fits 1.5 tube (closed)	2,000	1.3	3.3
120011	Standard	Semi Point	Single Tip		1,600	2.5	6.4
120012	Standard	Needle Point					

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com