

Totemax® Blood Collection Tray

Keep off the floor, work top or bed

Unique phlebotomy tray design and patented hook feature allows the tray to be on or over most bed rails, foot-boards and handrails found in a typical patient room

- Made from food grade HDPE resins for a durable non-breakable design
- Deep tray design helps reduce spills



Hooks Lock
in place when not in use



Hangs on bed rails,
footboards and handrails

Item No.	Color		L x W x H in	L x W x H cm
120021	Blue	■	17 x 11 x 11	43.2 x 27.9 x 27.9

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com