



## GENMark DNA Ladders

### 50bp DNA ladder




**GEN-DM50**

### 100bp DNA ladder



**GEN-DM100**

### 1kb DNA ladder




**GEN-DM1000**

**Storage** :-20°C for 24 months

**Recommended Loading Volume** :5µl/well

Features	GEN-DM50	GEN-DM100	GEN-DM1000
Size Range	50-1500bp	100-1500bp	250-10000bp
Band Number	17	11	13
Tracking Dye	Orange G	Orange G and XyleneCyanolFF	Bromophenol Blue and XyleneCyanolFF



**Ready to Use DNA ladder**

Pre-Mixed with loading dye for direct gel loading

**Sharp bands**

**Quick reference-** enhanced bands

**Stable-** Room temperature storage for 6month

## GENAXY Range of DNA ladders

Catalogue No	Description	Pack Size
GEN-DM100	<b>GEN-MARK</b> 100bp DNA Ladder 500ul/vial	1 Vial
GEN-DM1000	<b>GEN-MARK</b> 1kb DNA Ladder 500ul/vial	1 Vial
GEN-DM50	<b>GEN-MARK</b> 50bp DNA Ladder 500ul/vial	1vial

**Contact us for more details**

## GENAXY SCIENTIFIC PRIVATE LIMITED

155-156, C Block, Community centre Plot No. 32 to 57, Behind Janak Cinema Janakpuri, New Delhi – 110058

Email: [technical@genaxy.com](mailto:technical@genaxy.com), Ph : +91-11-47619999 (100 Lines),