



CONTAINER SEAL

Genaxy Container Seal is a sealing tape designed for the vast range of bottles and containers used in laboratory applications. Its unique elastomeric characteristics allow it to stretch and conform around the irregular perimeters of containers lids and covers to form a tight seal. It will block water vapor loss, reduce the rate of gas exchange and prevent the entry of contaminants.

Cat. No.	Color	Tensile	Adhesion	Operating Temp.	Application Temp	Material Thickness	Total Thickness
GEN-CS-LW	White	10lbs/in width	28 oz/in width	-40°C - 77°C	10°C-77°C	4.5 mil	5.5 mil
GEN-CS-XLW	White	10lbs/in width	28 oz/in width	-40°C - 77°C	10°C-77°C	4.5 mil	5.5 mil

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com