



## MARKER PEN

**SOLVENT RESISTANT** : Resistant to most solvents including xylene, alcohol, and DMSO. May be used in the presence of acetone showing moderate resistance. Ideal for histology cassettes, microscope slides and cryogenic vials. Markings will not run on frosted/painted labware. Markings are smear-proof, odorless and fade-resistant.

**TWIN TIP**: This Sharpie<sup>®</sup> pen has a fine tip at one end and an ultra fine tip at the other end. It will write on most surfaces.

**TEENY WRITER**: A retractable Sharpie<sup>®</sup> pen with an extra fine tip. Recommended for writing on all our labels for samples that will be stored at room temperature and not exposed to chemicals and/or caustic agents.

**SHARPIE FINE TIP & EXTRA FINE TIP SHARPIE**: Sharp fine and extra fine tip for minute detailings.

Cat. Number

Description

Qty. / Pack

GEN-SRP-BK

Solvent Resistant Pens, Black.

1

GEN-TTIP-BK

Twin-Tip teeny writer , Black

1

GEN-TWTR-BK

Teeny Writer, Black

1

GEN-FTP-BK

Fine Tip Marker, Black

1

GEN-EFTP-BK

Extra Fine Tip Marker, Black

1

**Genaxy Scientific Pvt. Ltd .**

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com