



PETRI SEAL

Genaxy Petri Seal stretches and conforms to form a tight seal and is temperature resistant, waterproof and chemically resistant. It resists abrasion and puncture and blocks vapor loss while reducing the rate of gas exchange. Use PetriSeal with microtiter plates, Petri dishes and culture plates to prevent to entry of contaminants. The seal remains flexible down to -30°C and can withstand continuous exposure from -40°C to +77°C.

Cat. No.	Color	Tensile	Adhesion	Operating Temp.	Application Temp	Material Thickness	Total Thickness
GEN-P-108-C	Clear	10lbs/in width	28 oz/in width	-40°C - 77°C	10°C-77°C	4.5 mil	5.5 mil
GEN-P-108-W	Clear	10lbs/in width	28 oz/in width	-40°C - 77°C	10°C-77°C	4.5 mil	5.5 mil

Genaxy Scientific Pvt. Ltd .