



TOUGH TAG

Genaxy Tough-Tags are chemically inert labels that strongly adhere to microtubes, test tubes, storage boxes and other microscope slides. They are temperature resistant from 121°C to -40°C and withstand autoclaving, boiling water baths, freezer temperatures, organic solvents and caustic agents. For best adhesion it is recommended that siliconized surfaces be avoided. ToughTags are perfect for accurately identifying important samples stored in tubes, micro plates or boxes that will be autoclaved or stored in moderate freezer temperatures

Cat. No.	Product Description	Pack Size
GEN-TTAG-SW	Tough -Tags White, 0.94" x 0.5" (L x W)	1000 tag
GEN-TTAG-SR	Tough -Tags Red, 0.94" x 0.5" (L x W)	1000 tag
GEN-TTAG-SB	Tough -Tags Blue, 0.94" x 0.5" (L x W)	1000 tag
GEN-TTAG-SG	Tough -Tags Green, 0.94" x 0.5" (L x W)	1000 tag
GEN-TTAG-SY	Tough -Tags Yellow, 0.94" x 0.5" (L x W)	1000 tag
GEN-TTAG-LW	Tough -Tags White, 1.28" x 0.5" (L x W)	1000 tag
GEN-TTAG-LR	Tough -Tags Red, 1.28" x 0.5" (L x W)	1000 tag
GEN-TTAG-LB	Tough -Tags Blue, 1.28" x 0.5" (L x W)	1000 tag
GEN-TTAG-LG	Tough -Tags Green, 1.28" x 0.5" (L x W)	1000 tag

Genaxy Scientific Pvt. Ltd .