



LABEL STATION

Label Station is used to organize labels for easy accessibility. Portable label station holds up to six standard size boxes of Tough-Tags[®], Cryo-Babies[®], and Tough-Spots[®]. Acrylic station includes a writing surface, a cutting edge for easy dispensing, and non-skid rubber feet. Counter-sunk holes provided for wall mounting. Acrylic station includes a writing surface, a cutting edge for easy dispensing, and non-skid rubber feet. Counter-sunk holes provided for wall mounting
Translucent blue 280 x 120 x 89 mm

Sr. No.	Catalog No.	Product Description	Pack Size
1	HS23460	Label Station, 280 x 120 x 89 (L x W x H)	1

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com