



CRYO BABIES AND CRYO TAGS

Cryogenic labels formatted for handwriting. Cryo-Babies[®] and Cryo-Tags[®] withstand both long term cryogenic storage conditions and conventional freezer storage. Cryogenic labels are chemically inert with a solvent resistant adhesive that will adhere to cryogenic vials, microtubes, test tubes, glass vials, beakers and metal racks. These can withstand flash/snap freezing in liquid nitrogen and withstand temperatures up to 80°C.

Cat. No.	Color	Size	Ideal Use	Temperature Range (°C)	Qty./Pack
GEN-CB-1000	White	1.28 x 0.50"	1.5ml tubes	-196 to 80	1000/roll
GEN-CT-2000	White	1.50 x 0.75"	General Purpose	-196 to 80	1000/roll
GEN-CB-0500	White	0.94 x 0.50	0.5ml tubes	-196 to 80	1000/roll

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com