



## WRITE PROTECT

A Clear Adhesive Overwrap to Protect Writing. It Defend writing against organic solvents and caustic agents. Comes in a handy dispenser and is perforated to fit 0.5ml and 1.5ml tubes. It can easily adhere to all plastics, glass and metals. Its temperature Ranges from -196°C to 150°C.

Sr. No.	Catalog No.	Dimension	Ideal Use	Temperature	Pack Size
1	GEN-WP-S	1.25" x 0.50" (L x W)	0.5 & 1.5ml	-196°C to 150°C.	1
2	GEN-WP-L	1.60" x 0.63" (L x W)	0.5 & 1.5ml	-196°C to 150°C.	1

**Genaxy Scientific Pvt. Ltd .**

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058  
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com