



GLACIERBRAND ICE BUCKET

Durable - Made of triple density expanded polystyrene creating a durable container that is virtually unbreakable

Super Insulating - Preserves ice up to 30% longer than conventional ice pans

Stackable - Designed with an inner edge enabling the pans to stack on top of one another for easy storage and to maximize bench-top space

Lightweight - Extremely lightweight (1.1oz, 4.7oz. and 11.4oz) and easy to handle. Feature a molded lip around the edge for carefree moving/carrying

Chemically Resistant - Resistant* to mild chemicals as well as ice-alcohol slurry, dry ice and liquid nitrogen (*Not recommended for long-term exposure to harsh chemicals or corrosives)

Sr. No.	Cat. No.	Bucket Size	Volume	Color	Pack Size
1	GEN-IB-2.5	10.25" x 4.75" (26.9cm x 11.1cm)	2.5 L	Blue	Each
2	GEN-IB-5.0	10.625" x 6.250" (26.9cm x 15.8cm)	5 L	Blue	Each

Genaxy Scientific Pvt. Ltd .