



GLACIERBRAND ICE PAN

Durable - Made of triple density expanded polystyrene creating a durable container that is virtually unbreakable

Super Insulating - Preserves ice up to 30% longer than conventional ice pans

Stackable - Designed with an inner edge enabling the pans to stack on top of one another for easy storage and to maximize bench-top space

Lightweight - Extremely lightweight (1.1oz, 4.7oz. and 11.4oz) and easy to handle. Feature a molded lip around the edge for carefree moving/carrying

Chemically Resistant - Resistant* to mild chemicals as well as ice-alcohol slurry, dry ice and liquid nitrogen (*Not recommended for long-term exposure to harsh chemicals or corrosives)

Sr. No.	Cat. No.	Pan Size	Volume	Color	Pack Size
1	GEN-IP-1	6.13" x 6.13" x 3.61" (15.6cm x 15.6cm x 9.17cm)	1L	Icy Blue	1
2	GEN-IP-4	12.25" x 8.75" x 4.75" (31.12cm x 22.2cm x 12.1cm)	4L	Icy Blue	1
3	GEN-IP-9	17.25" x 13.50" x 4.96" (43.8cm x 34.3cm x 12.63cm)	9L	Icy Blue	1

Genaxy Scientific Pvt. Ltd .