




Slide Mailers

Safe and secure transport of standard microscope slides
(76 x 26 mm and 3 x 1 in)




Cardboard Slide Mailers

Item No.	Color	Capacity	UOM	L x W x H in	L x W x H mm
HS9903	Natural 	1-Place	Pack of 25	3.7 x 1.7 x 0.2	95 x 42 x 5
HS9904	Natural 	2-Place	Pack of 36	4.1 x 3.1 x 0.2	103 x 79 x 5
120557	White 	4-Place	Pack of 20	4.7 x 5.7 x 0.3	120 x 140 x 6

- Heavy cardboard mailers for safe transport
- Cover that opens for viewing and thumb grooves for easy slide removal
- Securely holds slides with or without cover glass
- Flat surface area is ideal for labeling or writing necessary information



Plastic Slide Mailers




Item No.	Color	Capacity	UOM	L x W x H in	L x W x H mm
120556	Natural 	1-Place	Pack of 50	3.3 x 1.6 x 0.2	83 x 41 x 5
HS15983B	Natural 				
HS15983A	Assorted* 	2-Place	Pack of 25	3.7 x 2.7 x 0.2	93 x 68 x 5

*Assorted colors (Blue, Purple, Red, Orange, Yellow)

- Polypropylene protectors ideal for secure shipping and storing of slides
- Large, easy-to-secure snaps, provide safe storage for valuable slides
- Variety of colors brings some fun to the lab and enables color coding



Plastic 5-Place Slide Mailers

Item No.	Color	Capacity	UOM	L x W x H in	L x W x H mm
HS15982	Natural 	5-Place	Pack of 25	3.2 x 0.7 x 1.2	81 x 19 x 31
HS15983G	Assorted* 	5-Place	Pack of 25	1.7 x 0.9 x 3.5	43 x 24 x 88
HS15986	Natural 				

*Assorted colors (Blue, Purple, Red, Orange, Yellow)

- Polypropylene containers ideal for secure shipping, storing or staining slides
- Polypropylene containers hold five standard slides (76 x 26 mm and 3 x 1 in)
- Inside wall features slotted channels to keep slides separated and to make removal easier
- Lid closes securely to keep slides safe



CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com