

Freezer Thermometer

Protect frozen and chilled samples from potential spoilage

Monitor temperatures inside your freezer or refrigerator with the added convenience of sound and flashing LED lights to alert out-of-range temperatures

- Suction cup attaches the probe inside freezer/refrigerator while the thermometer secures to the outside with hook-and loop tape
- Large digital display shows freezer/refrigerator temperature plus ambient temperature
- Switches toggle for freezer/refrigerator monitoring and °F/°C
- **Set freezer range:** 14.9° to -9.4°F (-9.5° to -23°C);
Set refrigerator range: 35.6° to 46.4°F (2° to 8°C)
- Includes one AAA battery



Item No.	L x W x H in	L x W x H cm
HS28630	4.1 x .79 x 0.9	1.4 x 1.9 x 2.8

Probe Cable measures 0.9 m, Probe measures 5 mm OD x 17 mm L

**CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD**

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com