

Lead Rings

Stabilize vessels and labware on the bench or in water baths

- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- C-shape format allows some manual size variation for a customized fit



Item No.	Shape	Color	Flask Sizes mL	ID in	ID cm	OD in	OD cm	Weight g
HS8881A	C	Blue ■	125–500	1.9	4.2	3	7.5	215
HS8882A	O			1.9	4.8	3	7.7	210
HS8881B	C	Blue ■	250–1000	2	5.1	3.5	9	465
HS8882B	O			2	5.1	3.5	9	490
HS8881C	C	Blue ■	500–2000	2.1	5.4	3.9	9.8	595
HS8882C	O			2.2	5.7	3.9	9.8	610
HS8881D	C	Blue ■	1000–4000	2.6	6.6	4.4	11.3	890
HS8882D	O			2.8	7	4.5	11.5	935

**Weights may vary by +/- 3%*

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com