

Aluminum Weighing Boats & Dishes

Dishes for general-purpose weighing, dispensing, storage, food processing, evaporating or drying

Aluminum Weighing Boats

- Sturdy yet lightweight boats feature a smooth, flat bottom for sample handling and a flanged, crimped wall for rigidity
- Finger tab handle can be etched for identification prior to weighing

Aluminum Weighing Dish

- Designed for moisture determination balances with slightly tapered walls for stacking



Item No.	Size	Dia. x H <i>in</i>	Dia. x H <i>mm</i>	Volume <i>ml</i>	Packaging
HS14521A	Small	1.7 x 0.5	43 x 12	20	100
HS14521B	Medium	2.2 x 0.6	57 x 14	39	100
HS14521C	Large	2.8 x 0.6	70 x 16	50	100
HS14522	Large	4 x 0.3	101 x 8	30	50

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com