

smartSpatula®

Designed to fit perfectly with standard
0.2, 0.5, 1.5/2.0, 15 and 50 mL tubes


- Made of virgin, food-grade polypropylene
- Ideal for dispensing or sample collection of powders or soft solids
- Will not shatter when used with liquid nitrogen frozen samples
- Hollow core minimizes heat transfer—good for hot and cold liquids
- Antistatic, sterile and Eco-Friendly versions available
- Complies with FDA standards for food contact

Lab tools that scoop, stir and scrape.



STERILE EO

Always there
and always
clean!



Item No.	Description	Color		Suggested Tube mL	Packaging	L x W in	L x W cm	Spoon (Flare) in	Spoon (Flare) mm	Scoop/Pick L x W in	Scoop/Pick L x W mm	Taper
a 120005	Micro Antistatic	Opaque	☐	0.2 and 0.5	300/Bulk	5.5 x .2	14 x .4	—	—	4.3	11 (V-Tip)	15°
b 120008	Eco-Friendly	Opaque	☐	15 and 1.5	300/Bulk	5.5 x .3	14 x .7	9.8 x 5.5	25 x 14	11.8 x 1.2	30 x 3	
c 120003	Standard	Blue	■		300/Bulk							
	Standard	Opaque	☐	15 and 1.5	300/Bulk							
d 120007	Std. EtO Sterile	Opaque	☐		100/10 x 10pk	8.3 x .3	21 x .7	9.8 x 5.5	25 x 14	11.8 x 1.2	30 x 3	
	Std. EtO Sterile*	Opaque	☐		100/10 x 10pk							
e 120006	Macro	Green	■	50	150/Bulk	12.2 x .4	31 x 1	16.5 x 9.3	42 x 23.5	11.8	30	35°

*Non Toxic/Non Pyrogenic smartSpatula™ is a trademark of LevGo Inc.

CORPORATE OFFICE GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com