

## Lab Drying Rack

Bringing color to a mundane task

Most Versatile Rack in the market

- Channels designed for more useful configurations
- Large and Small pegs for all glassware varieties
- One of a kind Drying Basket to collect smaller items

Designed to protect against breakage

- Peg angle allows for better hold - especially large, heavy pieces
- Patent Pending SureGrip Hook securely grips glassware into place

Included:

- 72 Pegs ~ 36 Small (Green) ~ 36 Large (Blue)
- 2 Drying Baskets
- Drainage Tubing
- Mounting Hardware




### Drying Basket

to keep small items and slotted to hold pipets safe while drying



**SureGrip Hook**  
Patent-pending design securely grips glassware into place



Use on counter or wall  
Stands alone on integrated hinged stand. Mounts easily on wall with included guide and hardware

Can hold up to 77 lbs  
(34.9 kg)

Rack made of durable polystyrene

### Lab Drying Rack

Item No.	L x W x H in	L x W x H cm
120791	17.4 x 26 x 5.5	44.1 x 66 x 13.9

### Replacement Parts

Item No.	Description
120793	Peg 10 Piece Kit, (5) Small, (5) Large
120792	Drying Basket Kit, 2 ea
120813	Replacement Draining Set, Tube and Elbow

**CORPORATE OFFICE**  
**GENAXY SCIENTIFIC PVT. LTD**

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com