

Workstation Lab Mat

Safely facilitate workflow processes

- Protect counter tops from hot items up to 200°C (392°F)
- Lip around the edge can hold approx. 250 mL of liquid
- Small footprint, 12 x 14.9 x 0.1 in (30.5 x 37.9 x 0.2 cm)



Gridded Surface
for tube staging, preventing movement of objects

Reversible

Side One: combination gridded and smooth surface
Side Two: smooth surface only



Item No.	Color	L x W x H in	L x W x H cm
120747	Blue	12 x 14.9 x 0.1	30.5 x 37.9 x 0.2
120748	Purple		

Silicone Lab Mat

Large working space

The lab mat can be used with items such as centrifuges, vortex mixers, hot plates, stirrers, pipettes and so much more

- Mat features designated magnetic area to keep stir bars in place
- Protect countertops from hot items up to 200°C (392°F)



Raised Edges
enable liquid retention of approximately 300 mL of fluid



Large Working Surface

accommodates benchtop equipment, weigh boats, stir bars, scales, etc.



Item No.	Side #1 Color	Side #2 Color	L x W x H in	L x W x H cm
120506	Yellow	Grey/Blue	14 x 23.5 x 0.1	35 x 60 x 0.2
120507	Purple	Purple/Grey		
120599	Blue	White/Blue		

CORPORATE OFFICE GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com